





AutoSyringeTM Patent Pending

Irrigation ... Simplified Comfort | Convenience | Control

SKU: [408000] PRICE: \$595.00

- Eliminate hand fatigue associated with standard syringe irrigation
- Easy to fill, quick-change 20mL canisters
- Safe, controlled flow of irrigating solutions with variable flow rate settings.



AutoSyringe™ will **SAVE** you **Money, Time, Hand Fatigue, and Waste.**

Complete Kit

Includes:

(1) AutoSyringe™ Handpiece(20) Canister Liners(100) Plastic Barrier Sleeves

(3) 20mL Canisters

(3) Autoclavable Sleeves

(1) Charging Base

Additional Accessories:

 $\hbox{(2) Additional Canisters with Color Coding Stickers} \\$

[408001] \$30.00

(2) Autoclavable Sleeves

[408002] \$32.00

(100) Plastic Barrier Sleeves

[408003] \$12.00



INNOVATIVE SOLUTIONS

About Us

Since its beginnings in 1997, Vista Dental has produced quality dental products which offer superior clinical results. With our "Mountains Above the Rest" commitment to quality, creative product development, and impassioned commitment of its employees, Vista has emerged as a leading producer of brand name merchandise. Today Vista Dental designs, manufactures, packages and distributes more than 300 cataloged items, many of which offer significant advantages and margin improvement opportunities to both our end users and dealer partners.

Product Research and Development

In addition to a great number of trademarks, the company owns dozens of patents and patent pending applications related to a wide variety of unique dental devices and innovative dental equipment. Our team of engineers, scientists and clinicians provide a multi-disciplinary approach to evaluate processes and to improve upon current techniques and methodologies. We use state-of-the-art equipment and techniques to improve upon existing materials and we continue to partner with universities and dental professionals in the development of new and innovative products.

Our Dedication to Customer Satisfaction

Working with Vista Dental, you become part of the winning process:

We treat our customers with an overriding commitment to fairness and customer service.

We are committed to offering real product value at reasonable prices compared to industry competitors.

We strive to not only meet, but exceed, industry standards for quality systems, product reliability, and customer service.

Our success is measured in terms of customer satisfaction and brand recognition.

Results are achieved by serving our customers through faster, smarter, more aggressive selling systems, and by being more creative than our competitors.

Innovative Solutions for the dental industry

Leading the way in:

- Endodontic Irrigation Solutions
- Piezo Ultrasonic Tips & Equipment
- Whitening: In-Office & Take-Home
- Composite Delivery Systems
- Tissue Management Products
- · Comparison Product line: Vista vs Ultradent®*











COMPARE AND SAVE WITH VISTA'S ULTRADENT®* COMPARISON LINE!



tradent* Pri		<u>Ultradent®* Product</u>	VISTA
SAVE BIG		Tips & Syringes	Tips & Syringes
25%	pg. 7	Black Mini®* Brush Tips	Sol-u-Flo™
10%	pg. 8	Black Mini®* Tips	Volu-Flo™
60%	pg. 7	Blue Micro®* and Black Micro®* Tips	Pre-Bent Tips
60%	pg. 7	Endo-Eze®* Tips	3/4" Appli-Vac™ Bendable Tips
64%	pg. 9,35	Endo-Eze®* Tips	1" Appli-Vac™ Irrigating Tips
25%	pg. 8	Inspiral®* Brush Tips	Spira-Flo™
29%	pg. 8, 14	Metal Dento-Infusor®* Tips	Stat-Flo™
34%	pg. 8, 34	NaviTip®* Tips	Flexi-Glide™ Utility Tip
29%	pg. 34	NaviTip®* Tips	Visco-Tip
85%	pg. 10	Ultradent®* Luer Lock Syringes	Luer Lock Syringes
8%	pg. 7	White Mac™ & White Mini™ Laser Tips	Blue-Flo™
		Tissue Management & Impressions	Tissue Management & Impressions
6%	pg. 14	Astringedent®* Solution	Quick-Stat™ Liquid
35%	pg. 14	Ultrapak®* and Ultrapak®* E Cord	Stat-Pak™
25%	pg. 14	Ultrapak®* Packers	Vista-Pak™ Packers
45%	pg. 13	ViscoStat®* Clear Gel	Quick-Stat™ FREE
42%	pg. 13	ViscoStat®* Gel	Quick-Stat™ FS
		Prevention / Hygiene	Prevention / Hygiene
24%	pg. 17	Opalescence®* Smilebrush Toothbrushes	Vista Spring™ Toothbrushes
		Whitening	Whitening
38%	pg. 22	OpalDam®* Liquid Dam	Liquid Dam
29%	pg. 22	Opalescence®* Bosst PF Whitening	Fluorescent™ Doctors Kit Take-Home Whitening System
31%	pg. 21	Opalescence®* Whitening	Fluorescent™ Express White In-Office Whitening System
		Preparation	Preparation
20%	pg. 25	Consepsis®* Solution	CHX Restorative
31%	pg. 25	Seek®* and Sable™* Seek®* Solutions	Caries Indicator
12%	pg. 25	Ultradent®* Micro Applicators	Micro-Applicator Brushes
37%	pg. 26	Ultradent®* Porcelain Etch	Porcelain Etch
45%	pg. 26	Ultra-Etch®* Gel	Best-Etch™
		Endodontics	Endodontics
59%	pg. 33	ChlorCid®*Solution	Sodium Hypochlorite
38%	pg. 32	Consepsis®* Solution	CHX-Plus™
73%	pg. 37	TriAway™* Adapter	EZ-Lock™
29%	pg. 35	UltraCal®* XS Solution	Vista-Cal™
16%	pg. 36	Ultradent®* Citric Acid	Citric Acid
26%	pg. 31	Ultradent®* EDTA	EDTA
31%	pg. 36, 39	Ultradent®* Luer Vacuum Adapter	Micro Aspirators
		Equipment	Equipment
17%	pg. 51, 52	! Ultrawave®* Ultrasonics	Piezo Pilot™ & Co-Pilot™
		General Merchandise	General Merchandise
9%	pg. 11	Syringe Caps	Syringe Caps
*5.1.5.01		™ Whita Mini™ Saak® Sahla™ ChlorCid® Consensis® TriAway™ MayiTin® Ongloscanca® OnglDam® Ultrayer	



MOUNTAINS ABOUT THE REST INNOVATIVE PRODUCTS, TRUSTED VALUE



Smear OFF™

This patent pending 2-in-1 mix effectively replaces two commonly used solutions: EDTA and CHX. Removing both smear layer and bacteria!

Page: 33



Brite-N-Lite™ **Oral Illumination System**

This clip-on fiber optic LED device is autoclavable and easy to use with any aspirator.

Page: 63



ENDO ULTRA

The first and only cordless ultrasonic device with the power to create acoustic streaming and cavitation within the endodontic space.



Therma-Flo™ **Composite Warming Kit**

The heating technology of Therma-Flo™ significantly increases the flow characteristics of composites.

Page: 62



This patent pending device makes endodontic irrigation easy, convenient and affordable. This cordless, compact, battery operated device provides a consistent, controlled flow of irrigants.

Page: 61



Spira-Flo™ Brush Tips

Nylon brush fiber for controlled delivery of viscous materials, etchants and resins, with a curved cannula for optimal entry.





Tips & Syringes	6	
Tissue Management & Impressions	12	
Prevention / Hygiene	16	
10/L-10	20	
Whitening	20	
Preparation	24	
Endodontics	28	
Evacuation	38	
Equipment	44	

General Merchandise

SHOP SMARTER COMPARE & SAVE!

A tip for every application... see our extensive offering on the following pages!



- Allows light pressure on sulculus
- Ideal for profound hemostasis
- Internal spiral channel
- Minimizes bubbles
- Curved cannula for optimal entry
- Adjustable fiber brush

Compare to: Ultradent Save TO 64%



Appli-Vac 3/4"Tips™ Appli-Vac™ 1"Tips Appli-Vac™ Pre-tipped Syringes Blue-Flo™ Tips **Colored Luer Lock Syringes** Flexi-Glide™ Utility Tip **High Pressue Syringes Luer Lock Syringes** Micro Flowthru® Brush Tips Micro-Evac™ Tip Micro-Flo™ Tips NiTi SuperFlex™ Tips **Perioflex™ Irrigator Bottle Pocket-Probe**™ **Pre-Bent Tips** Safe-T™ Syringes Secure-Lock™ Appli-Vac™ Secure-Lock™ Vista-Probe™ **Sol-u-Flo™ Applicator Tips** Spira-Flo™ Brush Tips Stat-Flo™ Tips **US-12 Utility Syringe** Vista-Probe™ Tips Vista-Probe™ Pre-tipped Syringes Volu-Flo™ Tips

Appli-Vac™

Compare to Ultradent®*: Endo-Eze®* Tips

3/4" Bendable Tips

Our bendable 3/4" Appli-Vac™ and Secure-Lock™ are the ideal application tips for sealants, irrigating, aspirating, air/water delivery, flowable composites, etchants, resins, primers and more! Ideal tips for application and aspiration.

- · Blunt end, dead-soft needle resists crimping
- Extended 3/4" length for pockets, fistulas and surgical sites
- · Universal luer lock style hub



COLOR	GAUGE	20 / CUP	100 / BAG	120/PK
	16ga - White	[315016]	[315116]	
	18ga - Pink	[315018]	[315118]	
	19ga - Brown	[315019]	[315119]	
	20ga - Yellow	[315020]	[315120]	
	22ga - Black	[315022]	[315122]	
	25ga - Blue	[315025]	[315125]	
		\$6.40	\$21.10	
	3/4" Bendable Combo Kit			[315510]
	(20 tips of each gauge)			\$24.99

Pre-Bent Tips

Compare to Ultradent®*: Blue Micro® Tips and Black Micro®

Ideal for use with etchants, resins, and flowable composites

- Opaque photo-sensitive hub
- Blunt end pre-bent needle
- Universal luer lock hub



COLOR	GAUGE	20 / CUP	100 / CUP
	16ga - White	[315306]	[315316]
	20ga - Opaque	[315330]	[315320]
	22ga - Opaque Black	[315332]	[315322]
	25ga - Blue	[315335]	[315325]
		\$6.35	\$13.45

Blue-Flo[™] Tips

Compare to Ultradent®*: White Mac™* Tips and White Mini™* Laser Tips

Ideal for irrigation, aeration, and impression materials.

- Autoclavable
- 23ga soft plastic tip
- · Luer lock hub

20 / CUP	100 / BAG
[316202]	[316222]
\$8.70	\$42.40





Secure-Lock™ Appli-Vac™

3/4" Bendable Tips for increased safety and ease of use

- Blunt end, dead-soft needle resists crimping
- Extended 3/4" length for pockets, fistulas and surgical sites
- Universal luer lock style hub with Secure-Lock™ threads that lock tip in place



COLOR	GAUGE	20 / CUP	80 / CUP
	19ga	[314019]	[314119]
	20ga	[314020]	[314120]
	22ga	[314022]	[314122]
	25ga	[314025]	[314125]
		\$7.60	\$21.10

1/2" Bendable Applicator Tips

For etchants, resins and flowable composites

Luer lock hub



Color	GAUGE	20 / CUP	100 / CUP
	22ga - Black	[315202]	[315222]
	22ga - Orange	[315205]	[315225]
		\$5.99	\$12.60

Sol-u-Flo™ Applicator Tips

Ideal for use with low viscosity solutions

• 19ga tip with adjustable brush fiber length

- · Curved cannula for optimal entry
- Opaque luer lock hub



20 / CUP	100 / BAG	120/PK
[312202]	[312200]	[312205]
\$15.40	\$50.40	\$187.70



Micro Flowthru® Brush Tips

- Bendable brush tips bends up to 90° without crimping
- Non-linting, non-absorbent plastic applicator brush tips
- · Luer lock hub

Color	GAUGE	20 / CUP	100 / CUP
	20ga	[320820]	[321820]
	22ga	[320822]	[321822]
	25ga	[320825]	[321825]
		\$10.20	\$42.25



Spira-Flo™ Brush Tips Compare to Ultradent**: Inspiral** Brush Tips

Nylon brush fiber for controlled delivery of viscous materials, etchants and resins.

- Curved cannula for optimal entry
- Internal spiral channel
- · Adjustable brush tip fiber length
- Minimizes bubbles
- · Luer lock hub



100 / CUP	500 / BAG
[321400]	[321450]
\$60.50	\$230.15



Volu-Flo™ Tips

Compare to Ultradent®*:Black Mini®* Tips

Large 18ga tip opening to allow easy extrusion of viscous materials

- Designed to deliver large volumes of material
- Photo sensitive opaque tip maintains the flow of light-cured materials
- · Luer lock hub





Micro-Flo™ Tips

Patented

Perfect tip for magnification enhanced procedures. Ideal for controlled microscopic delivery of etchants and composites in most micro-restorative procedures.

- Disposable, polyimide tip provides unsurpassed delivery
- · Bendable needle offers areater access
- · Luer lock hub

COLOR	GAUGE	20 / CUP	100 / CUP
	Yellow	[406126]	[406626]
	Black	[406127]	[406627]
	Blue	[406130]	[406630]
		\$26.75	\$65.25



Stat-Flo[™] Tips
Compare to Ultradent*: Metal Dento-Infusor* Tips

Pre-bent 19ga Tips

Ideal for application of hemostatic agents

- Curved cannula for optimal visibility
- Padded-end brush tip




20 / CUP	100 / BAG	500/PK
[312102]	[312100]	[312105]
\$15.55	\$59.55	\$219.80



NiTi SuperFlex[™] Tips

SuperFlex[™] tips are slotted and side-vented for safe irrigation

- 2.5 X more flexible
- · Longer life cycle with less clogging
- 30ga nickel-titanium irrigating needle tips
- Adjustable needle angle for maximum flexibility
- AUTOCLAVABLE for cost-effective repeated use





Flexi-Glide™ Utility Compare to Ultradent*: NaviTip**

Patented

Our best choice for irrigation

- · Luer Lock Hub
- 30qa disposable (single use) tip
- Flexible, bendable needle facilitates access
- Tip can be cut to desired length



20 / CUP	50 / CUP
[406320]	[406350]
\$34.25	\$56.90



Appli-Vac™ Irrigating Tips Compare to Ultradent**: Endo-Eze** Tips

1" Bendable Irrigating Tips

Appli-Vac™ tips come individually

capped in 100/bag

- Slotted and side-vented for safe irrigation
- · Luer Lock Hub



2.02		456	
COLOR	GAUGE	20 / CUP	100 / BAG
	23ga	[315023]	[315523]
	27ga	[315027]	[315527]
	30ga	[315030]	[315530]
		¢7.05	624.00

^{*}Promotion is good only on 100/bag quantity.

Secure-Lock[™] Appli-Vac[™] Compare to Ultradent**: Lok-Tite** Feature

1" Bendable Irrigating Tips Luer lock hub with Secure-Lock™ threads to lock tip in place for increased safety and ease of use.

· Slotted and side-vented for safe irrigation



Appli-Vac[™] tips come individually capped in 80/bag







GET 1

Vista-Probe[™] Irrigating Tips

1" Bendable Irrigating Tips

capped in 100/bag




COLOR	GAUGE	20 / CUP	100 / BAG
	23ga	[315123]	[315623]
	27ga	[315127]	[315627]
	30ga	[315130]	[315630]
		\$21.55	\$90.50

^{*}Promotion is good only on 100/bag quantity.

Secure-Lock™ Vista-Probe™

1" Bendable Irrigating Tips Luer lock hub with Secure-Lock™ threads to lock tip in place for increased safety and ease

• Closed-end tip for side-port delivery



capped in 80/bag





COLOR GAUGE 20 / CUP 80 / BAG [314123] [314623] 23ga [314627] 27ga [314127] 30ga [314130] [314630] \$25.25 \$90.50

^{*}Promotion is good only on 80/bag quantity.



COLOR GAUGE 23ga [314023] [314523] 27ga [314027] [314527] 30ga [314030] [314530]

^{*}Promotion is good only on 80/bag quantity.

Appli-Vac™ Pre-tipped Syringes



Tipped with Appli-Vac™ **Irrigating Tips**

Features:

- Side-vented needles similar to Monoject® style
- Increase efficiency with our pre-tipped syringes.
- Pre-tipped syringes are ideal for easy irrigation of pockets, fistulas and canals.

GAUGE	QTY	3CC SYRINGE	12CC SYRINGE
23ga	100 / BOX	[316133]	[316113]
27ga	100 / BOX	[316137]	[316117]
30ga	100 / BOX	[316139]	[316119]
		\$42.05	\$46.30

Vista-Probe™ Pre-tipped Syringes



Features:

- Closed end, side-port needle.
- Increase efficiency with our pre-tipped syringes.
- Pre-tipped syringes are ideal for easy irrigation of pockets, fistulas and canals.



GA	JGE	QTY	3CC SYRINGE	12CC SYRINGE
23	ga	100 / BOX	[316233]	[316213]
27	ga	100 / BOX	[316237]	[316217]
30	ga	100 / BOX	[316239]	[316219]
			\$89.30	\$93.45

Luer Lock Syringes

Ideal for surgical site preparation.

- Textured finger grips to minimize slippage
- Graduated (cc) increments
- · AUTOCLAVABLE, non-sterile
- Standard luer lock style



Use with any Vista tips

QTY	1.2CC SYRINGE	3CC SYRINGE	6CC SYRINGE	12CC SYRINGE
100 / Pack	[316002]	[316030]	[316060]	[316012]
	\$28.70	\$12.30	\$17.99	\$22.50

High Pressue Syringes



(Shown with syringe cap) Winged 12cc High-Pressure syringe is made of rigid, clear polycarbonate.

QTY	1.2CC SYRINGE	PRICE	
60 / Pack	[316525]	\$36.80	

Colored Luer Lock Syringes

A fast, easy way to organize and identify syringes.

- Increases safety reduces incidence of syringe swap
- Available in 3cc and 12cc with graduated increments
- AUTOCLAVABLE and non-sterile





Color (80/Box)	QTY	3CC SYRINGE	12CC SYRINGE
	80 / BOX	[317158]	[317058]
	80 / BOX	[317159]	[317059]
	80 / BOX	[317160]	[317060]
	80 / BOX	[317161]	[317061]
		\$11.85	\$21.65
Assorted Pack		[317063]	[317062]
(20 per color)		\$10.65	\$19.45

Micro-Evac™ Tip


Patented

Ideal for removing moisture from canals.

- 25ga flexible, narrow tip
- Easily navigates curved canals
- Reduces the need for paper points
- Secure-Lock™ hub for greater security
- Includes (1) FREE Luer HVE Adapter



20 / CUP	50 / CUP
[406720]	[406750]
\$23.95	\$45.25



Syringe Caps

Compare to Ultradent®*: Syringe Caps

For preserving and sealing syringes and materials.

- White plastic caps with easy twist-open wing tabs
- Fits standard luer lock syringes



	25 / BAG	PRICE
Syringe Caps	[316910]	\$6.85

Pocket-Probe™

Patented

Ideal tips for irrigation of inter-dental spaces and perio pockets.

- 21ga tip
- Multiple perforations for thorough flushing
- Closed end, flexible 2-piece bent plastic irrigation tip

100 / CHB

• Universal luer lock style hub

20 / CHB

207 COI	100 / COI
[316120]	[316102]
\$27.05	\$83.30





Flexible 2-piece design

Use with Perioflex™ Irrigator Bottle. Perfect for patient take-home care.

Perioflex™ Irrigator Bottle





SKU	
[316200]	

US-12 Utility Syringe



Safe-T™ Syringes

Controlled, metered delivery

Ideal syringe for greater accuracy and control.

- · Ratcheted plunger emits a tactile and audible signal for each 0.2mL delivered
- · Easy, one-hand use
- AUTOCLAVABLE
- Universal luer lock

QTY	6CC SYRINGE	12CC SYRINGE
50 / Box	[316470]	[316480]
	\$33.70	\$39.35





Quick-Stat™ FREE

Quick-Stat™ TM Bulk Basic Kit
Quick-Stat™ FS

Quick-Stat™ TM Bulk Basic Kit
Quick-Stat™ TM COMPLETE Kit
Impression Trays
Vista-Pak™ Packers
Quick-Stat™ Liquid
Stat-Flo™ Tips
Impression Syringe
Stat-Pak™
Stat-Pak™+E

TISSUE MANAGEMENT & IMPRESSIONS

Quick-Stat™ FREE

25% Aluminum Chloride Clear Hemostatic Gel

- Surfactant lowers surface tension for better wound penetration.
- · Controls bleeding fast WILL NOT stain tissue
- Viscous gel will not flow onto surrounding tissue
- Includes Stat-Flo™ applicator tips





	SKU	PRICE:
tandard Kit	[504600]	\$18.00
1) 1.2mL prefilled syringes & (8) Stat-Flo™ tips		
ulk Syringe Kit	[504650]	\$41.75
) 30mL prefilled syringe, (20) Stat-Flo™ tips,		
0) 1.2mL empty syringes, and (1) docking port		
ılk Refill	[504660]	\$29.90
) 30mL prefilled syringe and (1) docking port		
alue Pack	[504640]	\$102.75
0) 1.2mL prefilled syringes & (50) Stat-Flo™ tips		

Quick-Stat™ TM Bulk Basic Kit

15% Ferric Sulfate - Bulk Hemostatic Gel Kits Includes:

(1) 30mL prefilled syringe

(20) 1.2mL empty syringes

(40) Stat-Flo™ applicator tips

(1) docking port

(2) Stat-Pak[™] cords - sizes 0, 1



SKU	PRICE
[503558]	\$82.25





Quick-Stat™ FS

15.5% Ferric Sulfate Hemostatic Gel

- Will not flow onto surrounding tissue
- Effective tissue management
- Quickly controls bleeding
- Includes Stat-Flo[™] applicator tips





	SKU	PRICE:
Standard Kit	[502804]	\$18.00
(4) 1.2mL prefilled syringes & (8) Stat-Flo™ tips		
Bulk Syringe Kit	[503550]	\$41.75
(1) 30mL prefilled syringe, (20) Stat-Flo™ tips,		
(20) 1.2mL empty syringes, and (1) docking port		
Bulk Refill	[503555]	\$29.90
(1) 30mL prefilled syringe & (1) docking port		
Value Pack	[502800]	\$102.75
(40) 1.2mL prefilled syringes & (50) Stat-Flo™ tips		

Quick-Stat™ TM Bulk Complete Kit Compare to Ultradent®*: ViscoStat®* TM IndiSpense®* Kit

15% Ferric Sulfate - Bulk Hemostatic Gel Kits Includes:

(1) 30mL prefilled syringe (20) 1.2mL empty syringes (20) Stat-Flo™ tips

DDICE

(1) docking port

(4) Stat-Pak™ Cords

(1) cord packer





PRICE	-
\$107.60	
	A CONTRACTOR OF THE PARTY OF TH
	The second second
	uick-star
STREET IN	ESPRINGE
TS.Ph feets Sulfide II	Canada Ca
	~
	DI STATE DE

TISSUE MANAGEMENT & IMPRESSIONS

Stat-Pak"

Compare to Ultradent®*: Ultrapak®* Cord

Knitted Retraction Cord

- Knitted Stat-Pak[™] cords are color coded for easy identification of cord size.
- Cords place easily with no fraying and no memory.
- Conform to sulcular anatomy with gentle outward force.
- Cords remain in place even when touched by a bur or diamond.
- Each bottle contains 100 in / 245 cm



COLOR	SIZE	SKU	PRICE
	000	[302611]	\$12.20
	00	[302612]	\$12.20
	0	[302613]	\$12.20
	1	[302614]	\$12.20
	2	[302615]	\$12.20







Stat-Pak™ +E

Compare to Ultradent*: Ultrapak* E Cord

Braided Retraction Cord

Cord is Impregnated with Epinephrine.

- Easy to place
- Cord will not fray and has no memory
- · Conforms to sulcular anatomy with gentle outward force
- Cords remain in place, even when touched by a bur or diamond
- Each bottle contains 72 in / 176 cm



COLOR	SIZE	SKU	PRICE
	00	[302633]	\$13.85
	0	[302632]	\$13.85
	1	[302631]	\$13.85

Vista-Pak™ Packers

Compare to Ultradent®*: Ultrapak®* Packers

Thin design - ideal for placing knitted retraction cords.

- Plain edge: to place and compress cord
- Serrated edge: for ease of use with knitted cord











VISTA-PAK™ PACKER	CIRCLET	STANDARD
Serrated	[302702]	[302704]
Plain	[302701]	[302703]
	\$28.50	\$28.50

Impression Syringe

Tip bends to 45° angle Seamless barrel construction eliminates splitting.

QTY	SKU	PRICE
50 / BOX	[313001]	\$15.00

Stat-Flo™ Tips

Compare to Ultradent®*: Metal Dento-Infusor®* Tips

Pre-bent 19ga Tips

Ideal for application of hemostatic agents.

- Curved cannula for optimal visibility
- · Padded-end brush tip



20 / CUP	100 / BAG	500/PK
[312102]	[312100]	[312105]
\$15.55	\$59.55	\$219.80



Quick-Stat[™] Liquid Compare to Ultradent**: Astringede

25% Aluminum Sulfate Solution Less irritating than aluminum chloride solutions.

- Higher pH reduces trauma to tissue
- Pleasant tasting mint flavor
- · No adverse reaction to impression materials



SIZE	SKU	PRICE
40 mL	[506306]	\$31.95

TISSUE MANAGEMENT & IMPRESSIONS


Impression Trays

Rigid acrylic with key-way perforations and side wall slots to ensure material retention.

- Superior retention
- Occlusal surface space minimizer
- Posterior rimlock minimizes secondary flow or gagging
- Finger rest ledges for easier removal from mouth
- Packaged: 12 / bag
- · Available in Clear and Green

#1 Large Upper [311001] #1 Large Upper [310001] #2 Large Lower [311002] #2 Large Lower [310002] #3 Medium Upper [311003] #3 Medium Upper [310003] #4 Medium Lower [311004] #4 Medium Lower [310004] #5 Small Upper [311005] #5 Small Upper [310005] #6 Small Lower [311006] #6 Small Lower [310006] #7 Upper Right/ Lower Left [311007] #7 Upper Right/ Lower Left [310007] #8 Upper Left / Lower Right [311008] #8 Upper Left / Lower Right [310008] #9 Anterior [311009] #9 Anterior [310009]	GREEN IMPRESSION TRAYS	SKU	CLEAR IMPRESSION TRAYS	SKU
#3 Medium Upper [311003] #3 Medium Upper [310003] #4 Medium Lower [311004] #4 Medium Lower [310004] #5 Small Upper [311005] #5 Small Upper [310005] #6 Small Lower [311006] #6 Small Lower [310006] #7 Upper Right/Lower Left [311007] #7 Upper Right/Lower Left [310007] #8 Upper Left / Lower Right [311008] #8 Upper Left / Lower Right [310008] #9 Anterior [311009] #9 Anterior [310009]	#1 Large Upper	[311001]	#1 Large Upper	[310001]
#4 Medium Lower [311004] #4 Medium Lower [310004] #5 Small Upper [311005] #5 Small Upper [310005] #6 Small Lower [311006] #6 Small Lower [310006] #7 Upper Right/ Lower Left [311007] #7 Upper Right/ Lower Left [310007] #8 Upper Left / Lower Right [311008] #8 Upper Left / Lower Right [310008] #9 Anterior [311009] #9 Anterior [310009]	#2 Large Lower	[311002]	#2 Large Lower	[310002]
#5 Small Upper [311005] #5 Small Upper [310005] #6 Small Lower [311006] #6 Small Lower [310006] #7 Upper Right/Lower Left [311007] #7 Upper Right/Lower Left [310007] #8 Upper Left / Lower Right [311008] #8 Upper Left / Lower Right [310008] #9 Anterior [311009] #9 Anterior [310009]	#3 Medium Upper	[311003]	#3 Medium Upper	[310003]
#6 Small Lower [311006] #6 Small Lower [310006] #7 Upper Right/ Lower Left [311007] #7 Upper Right/ Lower Left [310007] #8 Upper Left / Lower Right [311008] #8 Upper Left / Lower Right [310008] #9 Anterior [311009] #9 Anterior [310009]	#4 Medium Lower	[311004]	#4 Medium Lower	[310004]
#7 Upper Right/ Lower Left [311007] #7 Upper Right/ Lower Left [310007] #8 Upper Left / Lower Right [311008] #8 Upper Left / Lower Right [310008] #9 Anterior [311009] #9 Anterior [310009]	#5 Small Upper	[311005]	#5 Small Upper	[310005]
#8 Upper Left / Lower Right [311008] #8 Upper Left / Lower Right [310008] #9 Anterior [310009]	#6 Small Lower	[311006]	#6 Small Lower	[310006]
#9 Anterior [311009] #9 Anterior [310009]	#7 Upper Right/ Lower Left	[311007]	#7 Upper Right/ Lower Left	[310007]
2,	#8 Upper Left / Lower Right	[311008]	#8 Upper Left / Lower Right	[310008]
\$9.55	#9 Anterior	[311009]	#9 Anterior	[310009]
		\$9.55		\$9.55





Both Clear and Green - as low as 77¢ ea.

Perfect for Same - Day Placement of Restoration

THE #1 CHALLENGE Make a GREAT 1st IMPRESSION



Your patients will thank you!

The key to making a GREAT 1st impression:

Vista's Quick-Stat™ FREE is for the quick control of bleeding necessary for quality restorations.

Tip #1 - Quality Restorations:

If you are using a self etch followed by a conventional hemostatic and are not meticulous in cleaning all of the preparation, it will be challenging to obtain quality adhesion, potentially resulting in leakage and staining underneath the restoration.

Solution:

Vista's **Quick-Stat™ FREE** is a non-staining 25% Aluminum Chloride clear hemostatic that will NOT jeopardize the quality of your restorations.



Vista products are available through your preferred dealer



Dual Arch Fluoride Trays

HVE Tongue Cleaners

Interproximal Brushes

Piezo Co-Pilot™

Piezo Pilot™

Pocket-Probe™

Single Arch Fluoride Trays

Tongue Cleaners

Vista-Breeze™ Toothbrushes

PREVENTION / HYGIENE



Vista Spring[™] Tooth Brushes

Adult and youth styles - available with custom imprinting.

- Easy-grip handles
- Soft DuPont®* bristles for gentle cleaning
- Oval angle head for easy, effective access in hard-to-reach areas
- · Comfortable, ergonomic rubber grips prevent slipping
- Six assorted colors: green, blue, red, purple, orange and light blue

Something to **SMILE** about . . .



Unique "flex-zone" handle relieves pressure on sensitive teeth and gums.

Buy 3 or More Boxes, Get FREE custom imprinting!

Vista-Spring[™] Toothbrushes

Compare to Ultradent®*: Opalescense®* Smilebrush Toothbrushes

Adult and youth styles - available with custom imprinting.



CUSTOM IMPRINTING
ADD \$7.20 PER BOX

٥r

BUY 3 OR MORE BOXES, GET FREE CUSTOM IMPRINTING!

Vista-Breeze[™] Toothbrushes

Five assorted colors: blue, red, green, purple, and pink (not shown).




	PKG	SKU	PRICE	AS LOW AS
Vista Spring™ Youth Size	40 / box	[302110]	\$23.95	.60¢
Vista Spring™ Adult Size	40 / box	[302106]	\$23.95	.60¢
Vista Breeze™ Youth Size	40 / box	[302104]	\$23.95	.60¢
Vista Breeze™ Adult Size	40 / box	[302102]	\$23.95	.60¢

Interproximal Brushes

To remove plaque and food particles in hard-to-reach interdental spaces

- · Soft bristles with strong wire to resist bending
- Individually wrapped for patient take-home convenience
- Ergonomic handle for patient comfort and control

COLOR	SIZE	(36 / PK)	PRICE
	XXS (0.7 - 0.8 mm)	[302137]	\$23.85
	XS (0.8 - 1.0 mm)	[302136]	\$23.85
	S (1.0 - 1.2 mm)	[302135]	\$23.85




PREVENTION / HYGIENE

Tongue Cleaners

Great promotional item when imprinted.

• Durable ridges enhance removal of surface plaque and odor causing bacteria



Individually packaged in assorted colors



72 / BOX PRICE [309001] \$75.05

BUY 3 OR MORE BOXES GET FREE CUSTOM IMPRINTING



HVE Tongue Cleaners

Ideal for use with high volume suction.

- Fits all high speed suction devices
- Safe for children and adults
- Helps prevent halitosis and root caries
- Remove and dispose of biofilm on the tongue
- Autoclavable

QTY	SKU	PRICE
25 / BAG	[403214]	\$14.00






Single Arch Fluoride Trays

Ideal for foam and gel fluoride solutions.

- Flexible trays adapt to the contour of all arches and teeth
- Extended locking grip tabs
- Color coded for identification



IZE (100 / BAG)	PRICE
nall	[205101]	\$16.95
dium	[205102]	\$16.95
irge	[205103]	\$16.95
	mall dium	dium [205102]


TA SPRING™ TOOTHBRUSHE

Dual Arch Fluoride Trays

Ideal for foam and gel fluoride solutions.

- Color coded for easy identification.
- Flexible hinge with extended locking grip tabs.

· Flexible trays adapt to the contour of



COLOR	SIZE	(50 / BAG)	PRICE
	PEDO	[205204]	\$19.40
	SMALL	[205201]	\$19.40
	MEDIUM	[205202]	\$19.40
	LARGE	[205203]	\$19.40

BUY 3 GET 1 FREE BOX OF VISTA SPRING™ TOOTHBRUSHES*

Pocket-Probe™

Patented

Ideal for perio pocket irrigation.

- · Multiple perforations for thorough flushing
- Closed end, flexible 2-piece bent plastic irrigation tip
- Universal luer lock style hub

20 / CUP	100 / CUP
[316120]	[316102]
\$27.50	\$84.55





Use with Perioflex™ Irrigator Bottle. Perfect for patient take-home care.



ORDER NOW AND SAVE \$200.00



Pilot™

An exceptional range of power for a wide variety of perio, hygiene and endodontic procedures. Accepts most piezo ultrasonic tips.



Piezo Co-Pilot™ Dual Reservoir Scaler Features: (2) 300mL reservoirs with easy-to-see liquid levels.









Includes 8 Ultrasonic Tips



Acteon Style [402701] EMS Style [402711]

with **LED:** \$1,695

Acteon Style [402800]



Includes 10 Ultrasonic Tips

without **LED:** \$2,459

Acteon Style [402901] EMS Style [402911]

with **LED:** \$2,599

Acteon Style [402801]

Vista Dental Products

www.vista-dental.com

toll free 1.877.418.4782



Fluorescent™ Dr. Kit
Fluorescent™ EXpress WHITE™
Fluorescent™ Refill Kit
In-Office Whitening Liquid Dam Kit
Liquid Dam Refill Kit



Doctor supervised take-home whitening

Clear, easy-to`-use viscous take-home whitening gel formulated in 16% or 22% carbamide peroxide PLUS fluoride to help lower sensitivity.

BETTER whitening power with **LESS SENSITIVITY!**



Fluorescent™ Doctor's Kit **Take-Home Whitening System**

Compare to Ultradent®*: Opalescence®* Doctor Kit

Comes in a convenient zippered cloth pouch with tie-wrap holder for syringes. **Kit includes:**

- (8) 1.2mL Prefilled Syringes with Tips (2) EVS Tray Sheets 5" x 5" (.40mm)
- (1) Tray Case

	MINT	CITRUS	PRICE
16% Carbamide Peroxide	[504300]	[504310]	\$32.15
22% Carbamide Peroxide	[504400]	[504410]	\$32.15



FREE Marketing Materials with order!



- A) Three-sided display counter top display
- B) Vertical two-sided counter top display
- C) Q & A brochure: "Frequently Asked Questions"

Fluorescent™ Refill Kit Compare to Ultradent**: Opalescence** Refill Kit

Includes:

(40) 1.2mL prefilled syringes with tips packaged in (10) 4-packs.



	MINT	CITRUS	PRICE	
16% Carbamide Peroxide	[504350]	[504320]	\$124.45	
22% Carbamide Peroxide	[504450]	[504420]	\$124.45	



WHITENING



Whiten now... IN OFFICE!

35% hydrogen peroxide, non-light activated! Works incredibly fast!

Fluorescent™ EXpress White In-Office Whitening System

Compare to Ultradent®*: Opalescence®* Boost PF Patient Kit

System includes:

- (1) Dual-barrel Whitening Syringe
- (2) 1.2mL Desensitizer Syringes
- (1) 2mL Liquid Dam Syringe
- (5) 20ga Pre-bent Tips
- (2) Micro Brush Applicators
- (5) 20ga Micro Flowthru® Tips
- (1) White Surgical Evac Tube
- (1) Cheek Retractor
- (1) Whitening Shade Guide



Scan code to LEARN MORE about Fluorescent[®] EXpress White



PRICE

It's easy to whiten teeth in ONE office visit with Fluorescent™ EXpress White!

- **▶** COMPLETE Single procedure kit
- ▶ READY TO USE NO MIXING NEEDED!
- ▶ WORKS FAST In just ONE office visit!

Fluorescent™ EXpress White is formulated for FAST whitening in just one office visit. With 35% hydrogen peroxide, EXpress White removes tough stains and discolorations and is recommended for whitening teeth prior to placement of veneers and crowns.



FREE with order of

Fluorescent™ EXpress White In-office Whitening

Marketing Materials

Vertical two-sided and triangle counter top display



In-Office EXpress Whitening Kit	[504290]	\$49.50
Liquid Dam Kit	[316050]	\$41.00
(4) 2 mL Liquid Dam syringes		
(10) 20ga Pre-Bent tips		
(10) Vol-u-Flo™ Tips		
Liquid Dam Refill Kit	[316051]	\$139.00
(20) 2mL Liquid Dam syringes		

BEFORE AND AFTER

Fluorescent™ Take-Home Whitening is a dentist supervised system which is effective for removing many types of stains and discolorations due to every day exposure to food, beverages and life style habits. It is also recommended for whitening prior to placement of veneers and crowns.



BEFOREWhitening with
Fluorescent™
EXpress White.



AFTER
Whitening in
just one office
visit.



Superior Restorations

Vista's complete line of restorative products set a new standard by providing advanced control and predictability while offering dentists reliable and affordable products.







Vista Dental combines research and innovation to deliver some of the best solutions in Restorative dentistry.

Compare to:
Ultradent®*
Save TO 45%



Vista Dental Products

www.vista-dental.com

toll free 1.877.418.4782



Best-Etch™
Caries Indicator
CHX Restorative
Micro-Applicator Brushes
Micro-Restorative Best-Etch™
Porcelain Etch

PREPARATION

Caries Indicator

Compare to Ultradent®*: Seek® and Sable™* Solution

Vista-Red™ & Vista-Green™ Glycol-based solution to facilitate removal of carious dentin. Ensures maximum amount of inner dentin is preserved and infected dentin is removed. Use with Sol-u-Flo™ applicator tips.



	COLOR	INCLUDES	SKU	PRICE
Standard Pack		(4) 1.2mL prefilled syringes & (20) Sol-u-Flo™ tips	[502700]	\$30.90
Standard Pack		(4) 1.2mL prefilled syringes & (20) Sol-u-Flo™ tips	[502750]	\$30.90
Refill Value Pack		(4) 1.2mL prefilled syringes	[502705]	\$22.30
Refill Value Pack		(4) 1.2mL prefilled syringes	[502755]	\$22.30
Value Pack		(20) 1.2mL prefilled syringes	[503450]	\$62.15
Value Pack		(20) 1.2mL prefilled syringes	[503500]	\$62.15
Unit Dose		(30) individually packaged vials & brush applicators	[317009]	\$40.80
Unit Dose		(30) individually packaged vials & brush applicators	[317019]	\$40.80
Caries 10 mL Kit		10mL bottle, (6) mixing wells, and (24) Micro-Applicator brushes	[503409]	\$37.90
Caries 10 mL Kit		10mL bottle, (6) mixing wells, and (24) Micro-Applicator brushes	[503419]	\$37.90
Mixing Wells		100 / BAG	[503414]	\$13.70







CHX Restorative

Compare to Ultradent®*: Consepsis®*

2% Chlorhexidine Gluconate - Use prior to restorative preparations of crowns and composites

- Aids post-operative sensitivity to pulpitis
- Contains no surface modifiers
- Helps increase bond strength



	360	PRICE
16oz. (480mL) - Bottle	[317010]	\$58.15
30mL - Bulk Syringe	[503750]	\$26.80
Restorative Kit (4) 1.2mL prefilled syringes & (40) Sol-u-Flo™ tips	[502450]	\$22.30







Micro-Applicator Brushes


Compare to Ultradent®*: Micro Applicators

For etchants, sealants, bonding agents and hemostatic agents.

- Lint-free brush tip
- Non-absorbent
- Bendable
- As Low as .10¢ each



COLOR	SIZE	QTY	SKU	PRICE
	Ultra-Fine	100 / BOX	[320805]	\$13.30
	Fine	100 / BOX	[320800]	\$10.20
	Regular	100 / BOX	[320801]	\$10.20
	Regular	100 / BOX	[320802]	\$10.20
	Regular	100 / BOX	[320803]	\$10.20
	Regular	100 / BOX	[320804]	\$10.20



Best-Etch[™] Compare to Ultradent°*: Ultra-Etch°* Gel

37% Phosphoric Acid Etch

Designed to etch and prepare dentin and enamel.

- Viscous solution that stays in place
- Rinses clean with water and leaves no residue
- Pre-bent applicator tips ensure precise delivery
- Includes 22ga Pre-Bent Tips (25qa Pre-Bent tips are available on request.)
- Available in unit dose.

	INCLUDES	SKU	PRICE
Standard Kit	(4) 1.2mL prefilled syringes & (25) 22ga pre-bent tips	[502115]	\$16.30
Bulk Etch Kit	(1) 30 mL prefilled syringes, (50) 22ga pre-bent tips,	[502235]	\$42.60
	(5) empty 6cc syringes, & (1) docking port		
Bulk Syringe	(1) 30 mL prefilled syringe & (1) docking port	[502230]	\$31.20
Value Pack	(12) 3 mL prefilled syringes & (75) 2ga pre-bent tips,	[502200]	\$35.30
Unit Dose	(30) 0.6mL prefilled syringes & (30) 22ga pre-bent tips	[502225]	\$31.20





Porcelain Etch

Compare to Ultradent®*: Porcelain Etch

Viscous 9.5% Hydrofluoric Acid Solution. For porcelain veneers and inlays prior to cementation or repairs.

- Viscous solution
- Rinses clean away
- Silane solution used after etchant enhances the bonding of resin material

	INCLUDES	SKU	PRICE
Standard Kit (2) 1.2mL prefilled etch syringes		[502500]	\$33.25
	(2) 1.2mL prefilled silane syringes		
	(25) Spira-Flo™ brush tips		
Porcelain Etch Refill Kit	(4) 1.2mL prefilled syringes	[502510]	\$35.45
	(25) Spira-Flo™ brush tips		
Silane Refill	(4) 1.2mL prefilled syringes	[502505]	\$15.95



Micro-Restorative Best-Etch™

37% phosphoric acid etch for controlled micro-delivery of etchants.

- Designed for ideal application and enhanced delivery during microscope enhanced procedures.
- Comes complete with patented 30ga Micro-Flo™ tip for microscope enhanced procedures.
- Disposable, polyimide tip with bendable needle for greater access and unsurpassed delivery.

			SKU	PRICE
Micro-Restorative Be	[503200]	\$52.50		
(30) 0.6mL prefilled syring	ges & (30) 30ga M	icro-Flo tips		
MICRO-FLO™TIPS	20 / CUP	PRICE	50/CUP	PRICE
26ga	[406126]	\$26.75	[406626]	\$65.25
27ga	[406127]	\$26.75	[406627]	\$65.25
30ga	[406130]	\$26.75	[406630]	\$65.25







WHY PAY MORE...

For canal cleansing solutions?

TESTED & PROVEN for Superior Results



CHX-Plus™

This is our patent pending, 2% chlorhexidine gluconate solution with powerful wetting agents and proprietary surface modifiers that improves the depth of penetration, reducing bacteria 10X faster than normal 2% chlorhexidine.

Page 32



EDTA

17% solution with a neutral pH, facilitates the removal of smear layer and dentin mud and prepares the dentin wall for better adhesion of filling materials.

Page 31



Chlor-XTRA™


Our patent pending 6% sodium hypochlorite formula is enhanced with powerful wetting agents and proprietary surface modifiers, allowing it to penetrate lateral canals and isthmuses by lowering the surface tension, yielding a much more powerful oxidizing agent.

Page 30



Compare the results:





SEMs at 2000x magnification at the same section depth, irrigated with SmearOff™ vs QMix** brand 2-in-1 irrigating solution. Note the preservation of peritubular cuffs with SmearOff™, demonstrating that while SmearOff™ completely removes smear layer it is gentler on dentin than the QMix** brand. ¹Independently confirmed by Nelson Labs; Time kill study protocol #STP0158.2

Smear OFF™

An EDTA based formula, enhanced with Chlorhexidine, which removes the smear layer and kills 99.9%** of bacteria. This mix eliminates steps and saves time on each procedure.

Page 33

** Independently confirmed by Nelson Labs: Time Kill Study protocol #STP0158.2





Appli-Vac™ 1"Tips

CanalClean™ AA & CA Kits

Chlor-XTRA™

CHX-Plus™

Citric Acid

Disposable MTA Carriers

EDTA

EndoSafe™

EZ-Lock™

EZ-Lock™ Adapter

Fill Station™

Flexi-Glide™ Utility Tip

Micro Aspirators

Micro-Evac™ Tip

NiTi SuperFlex™ Tips

Pre-Tipped EZ-Lock™ 3/4" Bendable

Precision Irrigator™

Rinse-n-Dry™

Secure-Lock™ Appli-Vac™

Secure-Lock™ Vista-Probe™

SmearOFF™

Sodium Hypochlorite

Syringe Warmer

Versa-Brush™

Visco-Tip™

Vista-Blue™

Vista-Cal™

Vista-Probe™ Tips



Research Based for Best Outcomes

Endodontic retreatment rates are up over 35%, in part, due to improper irrigation techniques and protocol. With advancements in file technology, practitioners are now passing fewer files into the canal than ever before. Irrigation is vital for reliable root canal procedures, but passing fewer instruments means doctors no longer use as much irrigant solution.¹

Source: http://www.aae.org/rootcanalspecialists/dentalprofessionalsandstudents/factsheet/

Everything you need, when you need it.

• 1 *Box* = 1 *Procedure*

CanalClean™ 2S Kit

for Complex Anatomy (multi-root canals)

Step 1 – 12cc Chlor-XTRA™ Enhanced Sodium Hypochlorite

Step② – 6cc SmearOFF™ 2-in-1 Solution EDTA enhanced with Chlorhexidine **Step**⑤ – 6cc Chlor-XTRA™ Enhanced Sodium Hypochlorite

+ (5) 30ga side-vented irrigating tips

+ (1) Micro Aspirator with tip

QTY	SKU	PRICE
4 kits	[504850]	\$32.00





See pg. 30

+ SmearOFF™ 2-in-1 Mix

See pg. 33

Enjoy the convenience of a two solution kit!

- ·New CanalClean™ 2S offers everything needed for optimal canal cleanliness, in a convenient 2 solution kit.
- ·CanalClean™ 2S will save you time, eliminating the need for EDTA, CHX, and a rinsing agent!
- ·Chlor-XTRA™ and SmearOFF™ are 100% compatible and WILL NOT form a precipitate when mixed, making irrigation simple & predictable.

CanalClean™ 4S Kit

for Anterior Anatomy

Step 1 – 12cc NaOCl – pharmacy controlled 6% concentration (Chlor-XTRA™)

Step 2 – 3cc EDTA 17% with optimized pH

Step 3 - 3cc Ethyl Alcohol Rinse

Step 4 – 3cc CHX-Plus™

Dry: Micro-Aspirator with tip included

QTY	SKU	PRICE
4 kits	[504800]	\$32.00





Ideal product for GP's doing occasional endo Simple. Convenient. Affordable.

Vista products are available through your preferred dealer

Chlor-XTRA™

Multiple Patents Pending

Vista's Chlor-XTRA™ is an enhanced 6% sodium hypochlorite solution designed for irrigation, debridement and cleansing of root canals during and after instrumentation. The addition of powerful wetting agents and proprietary surface modifiers enables it to penetrate into lateral canals and isthmuses by lowering the surface tension.



Clinically Proven:

- To work 4X FASTER than ChlorCid®*
- To reduce post-treatment complications
- To improve contact and dissolve more tissue
- To increase lubricity
- To decrease surface tension

Chlor-XTRA™ offers:

- Greater oxidizing power
- 2+ times more wetting ability
- 2+ times more digestive
- · Half the surface tension



Scan to see the comparison to standard Sodium Hypochlorite

BUY 4 GET 1 160Z. BOTTLE FREE* BUY 8 GET 1 FILL STATION™ FREE* (ONLY APPLIES TO 160Z. BOTTLES)

	SKU	PRICE
16oz. (480mL) bottle	[503800]	\$61.90
3cc Pre-filled syringes (individually packaged)	[503825]	\$15.25
12cc Pre-filled syringes (individually packaged)	[503850]	\$15.25
Kits include 10 pre-filled syringes and 27ag side-vent irrigating tipe	5	



Effectiveness of pulp and tissue dissolution in 30ml solution for 45 minutes





Start of the test: Equal tissue sample is based on weight (0.125g / 125mg). Each sample was in 30ml of solution (for 4.166mg of tissue per ml of solution).

Tissue dissolution by sodium hypochlorite: effect of concentration, temperature, agitation and surfactant

Aim: Sodium hypochlorite is the most commonly used endodontic irrigant due to its antimicrobial and tissue dissolving activity. The aim of this study was to evaluate and compare the effects of concentration, temperature and agitation on the tissue dissolving ability of sodium hypochlorite. In addition, a hypochlorite product with added surface active agent was compared with conventional hypochlorite solutions.

Methods: Three sodium hypochlorite solutions from two different manufacturers in concentrations of 1%, 2%, 4% and 5.8% were tested at room temperature, 37°C and 45°C with and without agitation by ultrasonic and sonic energy and pipetting. Distilled and sterilized tap water were used as controls. Pieces of bovine muscle tissue (68±3) were placed in 10 ml of each solution for five minutes. In selected samples, agitation was performed for one, two or four 15 sec periods per each minute. The tissue specimens were weighed before and after treatment, and the percentage of weight loss was calculated.

Results: Weight loss (dissolution) of the tissue increased almost linearly with the concentration of sodium hypochlorite. Higher temperatures and agitation considerably enhanced the efficacy of sodium hypochlorite. The effect of agitation on tissue dissolution was greater than that of temperature, continuous agitation resulting in fastest tissue dissolution. Hypochlorite with added surface active agent was most effective in tissue dissolution in all experimental situations.

J Endod. 2010 Sept; 36(9); 1558-1562 S Stojicic, DDS, MSc ,*† S Zivkovic, DDS, PhD,† W Qian, DDS, PhD*and M Haapasalo, DDS, PhD* *Division of Endodontics, Department of Oral Biological & Medical Sciences, Faculty of Dentistry, University of British Columbia, Vancouver, Canada; † Department for Restorative Dentistry and Endodontics, Faculty of Dentistry, University of Belgrade, Belgrade, Serbia.

Antimicrobial and Tissue Dissolving Capabilities of **Different NaOCI Concentration**

The objective of this in vitro study was to establish the antimicrobial efficacy and the potential to dissolve vital pulp using different concentrations of sodium hypochlorite (NaOCI) solutions. The following solutions were tested: (1) 2% NaOCI (Thick Jik); (2) 3.5% NaOCI (Jik); 3) 6% NaOCI (Vista); (4) 6% NaOCI (Dent-Ionics); (4) 6% NaOCI (Chlor-Xtra[™]) and (5) Water (Control).

No significant inhibition of E. faecalis was observed with sterile water (control) and 3.5% NaOCI. However, the 2% NaOCI and all the other 6% NaOCI solutions showed significant inhibition against E. faecalis. Chlor-Xtra™ dissolved the pulp tissue in significantly less time compared to all the other groups (p<0.05).

Chlor-Xtra™ offered the best antimicrobial properties and tissue dissolving capabilities.

P. VAN DER VYVER, F. BOTHA, and F.A. DE WET University of Pretoria, School of Dentistry, Pretoria, South Africa



COMPARATIVE ANALYSIS Chlor-XTRA™ ChlorCid®*			
Free Chlorine (%)	• 6	2.3	
Surface Tension (dynes / cm)	• 29	45.8	
Dissolution Rate (mg / min)	2.3	0.43	
Cost	60.12/ml	\$0.37/ml	

After 45 minutes Chlor-XTRA™ completely dissolved the tissue. In the ChlorCide* container, the sample tissue will not be completely dissolved for another THREE hours

EDTA

Compare to Ultradent®*: EDTA

17% Solution - Specifically engineered with LOWER pH for optimal performance! For best outcomes, a LOWER pH is more important than concentration. As pH increases, the availability of calcium ions from hydroxyapatite for chelation decreases . . . a pilot study using concentrations of alkaline (pH 11.3) solutions of EDTA showed them to be less effective in removing the smear layer.

Vista's EDTA:

- Prepares dentin walls for better ADHESION of filling materials
- Opens dentin tubules
- Removes smear layer and dentin mud
- Enhances cleansing





Scan to see the optimal pH for EDTA solutions should be between 7.5 - 9.

BUY 1 GET 1 160Z. BOTTLE CHLOR-XTRA™ FREE*
BUY 8 GET 1 FILL STATION™ FREE*

(ONLY APPLIES TO 160Z. BOTTLES)

	SKU	PRICE
30mL Bulk Syringe	[503700]	\$9.60
16oz. (480mL) bottle	[317011]	\$54.25
4oz. (480mL) bottle	[317001]	\$16.80
Unit Dose (20) 0.6mL Individually packaged prefilled syringes & (20) 27ga side-vent irrigating tips	[317002]	\$23.15

Rinse-n-Dry™

95% Ethyl Alcohol Solution
Penetrates lateral canals and isthmuses for improved drying and cleansing.

- Low surface tension
- · Use after smear layer removal with EDTA
- Clears canal of all sodium hypochlorite before the final irrigation step with CHX-Plus™

	SKU	PRICE
4oz. bottle	[317003]	\$9.65



Demineralizing effects of EDTA at different concentrations and pH.

The purpose of this study was to compare effects of concentration and pH variations of EDTA on dentin demineralization. Twenty extracted, human permanent teeth with single canals were used in this study. Demineralizing effects of EDTA solutions at 10% and 17% concentrations at pH 7.5 and 9.0 were determined by measuring the amount of liberated phosphorus 1, 3, 5, 10, and 15 min after exposure.

Results showed that the amount of phosphorus liberated from dentin was greater with increased EDTA concentration and increased time of exposure, and was more effective at neutral pH than pH 9.0. The pH of EDTA solutions did not display any significant alterations during the demineralization process.

J Endod. 2009 Jan;35(1):95-7. Epub 2008 Nov 7. Williamson AE, Cardon JW, Drake DR. Department of Endodontics, University of Iowa College of Dentistry, Iowa City, Iowa, USA. anne-williamson@uiowa.edu

A scanning electron microscope evaluation of four root canal irrigation regimens.

A scanning electron microscope was used to evaluate the debridement capabilities of four irrigation regimens on both instrumented and uninstrumented root canal surfaces.

The combination of NaOCI and EDTA used alternately completely removed the smear layer from the instrumented root canal surfaces as well as the pulpal remnants and predentin from the uninstrumented surfaces. In addition, the combination of NaOCI and EDTA caused the exposed calcospherites on the uninstrumented surfaces to have an eroded appearance.


J Endod. 1987 Apr;13(4):147-57. Baumgartner JC, Mader CL.


CAUTION

Not all EDTA Solutions are as effective as they claim!

SMEAR LAYER REMOVAL

SEM Images





Vista - 17% EDTA

11.56 pH - 18% EDTA

Studies Show:

Lower pH is CRITICAL to performance and **MORE** important than concentration*

50 years of research shows that pH matters!				
Compare: Ultradent®**				
pH Level:	8.5	11.56		
Concentration:	17.0%	18.0%		

pH Measurements: Using Fischer Scientific AR 20 model, Electrode standardized using pH Standard solutions 10 and 4 "Study data: Nikiforuk G, Sreenby L. Demineralization of hard tissues by organic chelating agents. Science 1951; 114:560 Cury JA, Bragotto C, Valdrighi L. The demineralizing efficiency of EDTA solutions on dentin. Hilleunce of pH. Oral Surg Oral Med Oral Pathol. 1981Oct; 52(4): 446-8 O'Connell MS, Morgan LA, Beeler WJ, Baumgartner JC. A comparative study of smear layer removal using different salts of EDTA J Endod. 2000 Dec; 26(12): 739-43 Serper A, Calt S. The demineralizing effects of EDTA at different concentrations and pH. J Endod. 2002 Jul; 28(7):501-2 Obst JJ. Chelation and dentistry: conceptual orientation. NY J Dent 1962; 30:331-5

Vista products are available through your preferred dealer

CHX-Plus™

Compare to Ultradentes: Consepsises Solutions

Multiple Patents Pending

Enhanced 2% Chlorhexidine Gluconate Solution Kills multispecies biofilm **2X faster** than standard 2% chlorhexidine. Contains 2% Chlorhexidine Gluconate plus powerful wetting agents and proprietary surface modifiers to improve the depth of penetration and kills bacteria 10X faster than normal 2% Chlorhexidine.



BUY 4 GET 1 160Z. BOTTLE FREE* BUY 8 GET 1 FILL STATION™ FREE* (ONLY APPLIES TO 160Z. BOTTLES)

	SKU	PRICE
16oz. (480mL) bottle	[503900]	\$58.15
Endo Kit	[502925]	\$13.70
(4) 1.2 mL prefilled syringes & (20) 27ga side-vent irrigating tips		
Endo Refill Pack	[503930]	\$20.70
(20) 1.2ml prefilled syringes		
30mL Bulk Syringe	[503650]	\$25.99
(1) 30mL prefilled bulk syringe & (1) docking port		
30mL Bulk Syringe Kit	[503400]	\$36.05
(1) 30mL prefilled bulk syringe, (20) 1.2cc syringes,		
(20) 27ga side-vent irrigating tips, and (1) docking port		
Prefilled Syringes Refill	[502455]	\$7.30
(4) 1.2mL prefilled syringes		
Unit Dose	[503925]	\$32.75
(30) 0.6mL Individually packaged prefilled syringes		
& (30) 27ga side-vent irrigating tips		


Evaluation of the effect of two chlorhexidine preparations on biofilm bacteria in vitro: a three dimensional quantitative analysis.

His study developed a biofilm model that closely mimicked in vivo biofilm to determine its susceptibility to endodontic antimicrobial irrigants by three-dimensional quantitative analysis. Collagen-coated hydroxyapatite (C-HA) and uncoated hydroxyapatite (HA) disks were inoculated with dispersed subgingival plaque for three weeks. Thick biofilms rich in spirochetes were formed on both substrates. Biofilms were subjected to 1-, 3-, and 10-minute exposures to CHX-Plus™ (Vista Dental Products, Racine, WI) and 2% chlorhexidine (CHX). After treatment, the volume ratio of dead bacteria to all bacteria in biofilms, indicated by the ratio of red and (red + green) fluorescence, was analyzed by confocal laser scanning microscopy for each medicament.

The proportion of killed bacteria depended on type of irrigant and time of exposure in both C-HA and HA biofilm models (p = 0.00). CHX-Plus™ showed higher levels of bactericidal activity at all exposure times than 2%CHX (p < 0.001). The C-HA biofilm was thicker than the HA biofilm. Less bacteria were killed in C-HA biofilm than in the HA model. J Endod. 2009 Jan;35(7):981-5. Epub 2008 Nov 7. Shen Y, Qian W, Chung C, Olsen I, Haapasalo M. Division of Endodontics, Department of Oral Biological and Medical Sciences, University of British Columbia, Vancouver, Canada.

Data reported on C-HA Biofilm

Antimicrobial effectiveness of Endo Solution (CHX-Plus™)



Test time: 3 minutes RED color represents the volume ratio of dead bacteria.

Study data:

Shen Y, Qian W, Chung C, Olsen I and Haapasalo M. Evaluation of the effect of two chlorhexidine preparations on biofilm bacteria in vitro: a threedimensional quantitative analysis. J Endod 2009; 35:981-5.

The synergistic antimicrobial effect by mechanical agitation and two chlorhexidine preparations on biofilm bacteria.

This study investigates whether mechanical agitation (ultra-sonic or sonic) improves effectiveness of chlorhexidine against biofilm bacteria in vitro. Collagen-coated hydroxyapatite (CHA) disks were exposed to dispersed subgingival plaque for 3 weeks at 37 degrees C. The multispecies biofilms established were subjected for 1 and 3 minutes to CHX-Plus™ (Vista Dental Products, Racine, WI) and 2% chlorhexidine (CHX), with or without mechanical agitation. After treatment, the amount of dead bacteria in biofilms was analyzed by viability staining and confocal laser scanning microscopy (CLSM). The morphology of biofilms, with or without mechanical agitation, was also examined by CLSM.The structure of the biofilm showed no obvious change when the surrounding solutions were exposed to continuous ultrasonic or sonic agitation. Combined use of mechanical agitation and chlorhexidine had a more pronounced anti-microbial effect against the biofilms than either one alone.

J Endod. 1987 Apr;13(4):147-57 Shen Y, Stojicic S, Qian W, Olsen I, Haapasalo M. Division of Endodontics, Department of Oral Biological and Medical Sciences, University of British Columbia, Vancouver, Canada.

Smear*OFF*™

Patent Pending

Proprietary EDTA / CHX mix.

Formulated with wetting agents and surface modifiers. For best outcomes, our patent pending 2-in-1 mix effectively replaces two commonly used solutions: EDTA and CHX.


The obvious choice over QMix®* brand


- Superior chelation
- Won't form a precipitate when mixed with NaOCl
- No additional rinse needed when used with NaOCI
- Better calcium suspension
- Save 33% over QMix** brand

SmearOFF™ removes both smear layer AND bacteria!

Results have shown that SmearOFF $^{\text{m}}$ is equally as effective in removing smear layer as QMix $^{\text{ex}}$ brand!

Compare the results:





SEMs at 2000x magnification at the same section depth, irrigated with SmearOff™ -vs-QMix®* brand 2-in-1 irrigating solution. Note the preservation of peritubular cuffs with SmearOff™, demonstrating that while SmearOff™ completely removes smear layer it is gentler on dentin than the QMix®* brand.

¹Independently confirmed by Nelson Labs;

Time kill study protocol #STP0158.2







	SKU	PRICE
16oz. (480mL) bottle	[317111]	\$92.85
4oz. (120mL) bottle	[317112]	\$30.30
Unit Dose	[317114]	\$39.00



Scan to see the scientific data on SmearOFF™ vs QMix®* Time Kill Study

Sodium Hypochlorite

Compare to Ultradent®*: ChlorCid®* Solution

3% and 6% Solutions

For the debridement and irrigation of root canals during and after instrumentation.

(10) prefilled syringes

(10) 27ga-vent irrigating tips



	3%	6%	PRICE
3mL syringes	[502355]	[502350]	\$10.75
12mL syringes	[502315]	[502300]	\$10.75
16oz. (480mL bottle)	[317018]	[317008]	\$12.85



Irrigation Promo Pack

A winning combination for optimal canal cleanliness and improved outcomes!

The only *TRUE* 2-step irrigation protocol

No rinse needed when using Chlor-XTRA™
in conjunction with SmearOFF™ 2-in-1!

Includes:

(1) Chlor-XTRA™ (16 oz. bottle)(1) SmearOFF™ (16 oz. bottle)

	SKU	PRICE
Irrigation Promo Pack	[317113]	\$79.00



Vista products are available through your preferred dealer

Flexi-Glide[™] Utility Tip Compare to Ultradent*: NaviTip** Tips Patented

Our best choice for irrigation.

- 30ga disposable (single use) tip
- Flexible, bendable needle facilitates access
- Tip can be cut to desired length
- · Luer Lock Hub



NiTi SuperFlex™ Tips

Patented

NiTi SuperFlex™ The "Gold Standard" for endodontic irrigation tips.

- 2.5 X more flexible
- · Longer life cycle with less clogging
- 30ga nickel-titanium irrigating needles
- Adjustable needle angle for maximum flexibility
- ${\bf \cdot} {\it AUTOCLAVABLE} \ for \ cost-effective \ repeated \ use$



EndoSafe™

Patent Pending

Apical Irrigation / Evacuation

Unique bellows feature allows you to move from one canal depth to another without modifying cannula. Ratcheting Safe-T™ syringe offers greater control in delivery of irrigants. Effective, Safe, and Non-Clogging! Our non-clogging EndoSafe™ is a simple negative pressure, device providing simultaneous irrigation and evacuation, comparable in performance to the Discus EndoVac®*! The collapsing bellows enhances irrigation, allowing users to move easily from one canal depth to the next with no need to modify the irrigating cannula. This eliminates the chance of clogging and saves valuable time during procedures.

- Simultaneous irrigation and aspiration
- Fluid flows from apical to coronal region to virtually eliminate clogging
- Ratcheting feature allows for a tactile and audible signal for every
 0.2mL delivered
- Easy one-hand use no assembly needed
- Disposable
- Includes flexible 30ga tip and 12cc Safe-T™ syringe
- Ideal for use with Vista endodontic solutions.

Working Length Delivery

- + Contact Time
- + Effective Irrigation
- = Improved Outcomes

QTY	SKU	PRICE
10 / BOX	[503980]	\$52.00

*According to recent studies showing no significant difference between the EndoSafe™ and the EndoVac® in removing debris from the canal.





Visco-Tip™

Compare to Ultradent*: NaviTip* Tips
Patented

The ONLY TIP specially designed for easy delivery of Calcium Hydroxide.

- Bendable, flexible 25qa tip
- Secure-Lock[™] hub for safe delivery



QTY	SKU	PRICE:
20 / CUP	[312302]	\$27.95
50 / CUP	[312300]	\$61.75



Appli-Vac™ Tips

Compare to Ultradent®*: Endo-Eze®* Tips

1" Bendable Irrigating Tips

- Slotted and side-vented for safe irrigation
- · Luer Lock Hub

COLOR	GAUGE	20 / CUP	100 / BAG
	27ga	[315027]	[315527]
	30ga	[315030]	[315530]
		\$7.85	\$34.00









Secure-Lock™ Appli-Vac™

Compare to Ultradent*: Lok-Tite** Feature

1" Bendable Irrigating Tips

Luer lock hub with Secure-Lock™ threads to lock tip in place for increased safety and ease of use.

• Slotted and side-vented for safe irrigation

COLOR	GAUGE	20 / CUP	100 / BAG
	27ga	[314027]	[314527]
	30ga	[314030]	[314530]
		\$9.45	\$33.90

Vista-Probe™ Tips

1" Bendable Irrigating Tips Closed-end tip for side-port delivery

• Luer Lock Hub









individually capped in 100 / bag



Appli-Vac[™] tips come individually capped in 100/bag

Secure-Lock™ Vista-Probe™ Compare to Ultradent®*: Lok-Tite®* Feature 1" Bendable Irrigating Tips Luer lock hub with Secure-Lock™ threads to

Luer lock hub with Secure-Lock™ threads to lock tip in place for increased safety and ease of use.

· Closed-end tip for side-port delivery

COLOR	GAUGE	20 / CUP	80 / BAG
	27ga	[314127]	[314627]
	30ga	[314130]	[314630]
		\$25.25	\$90.50

Micro-Evac™ Tip

Compare to Ultradent®*: Capillary Tip

Patented

- 1" Bendable Irrigating Tips Ideal for removing moisture from canals.
- Slotted and side-vented for safe irrigation
- Luer lock hub with Secure-Lock™ threads lock tip in place for increased safety and ease of use.



QTY

20 / CUP

50 / CUP





Ideal replacement for **Ultradent®*** Capillary Tip

Disposable MTA Carriers

- For one-time use
- Delivers MTA directly to the apex of the canal
- · Bendable cannula
- Autoclavable

SIZE	QTY	SKU	PRICE
0.80mm	5 / PACK	[807017]	\$51.25
0.99mm	5 / PACK	[807016]	\$51.25
1.60mm	5 / PACK	[807015]	\$51.25
Variety Pack	5 / PACK	[807019]	\$51.25
(2) 1.6mm, (1) 0.80mm, (2) 0.99mm			



BUY 2 GET 1 FREE VARIETY PACK

Vista-Cal™

Compare to Ultradent®*: UltraCal®* XS Paste

PRICE

\$23.95

\$45.25

12.5 pH aqueous solution

SKU

[406720]

[406750]

Ideal for temporary filling and cleansing root canals.

- Contains no more than 35% Calcium Hydroxide paste
- Includes patented Visco-Tip™ and high pressure syringes for the most effective delivery of Calcium Hydroxide.



	SKU	PRICE
Vista-Cal™ Standard Kit	[502600]	\$38.75
(10) Visco-Tips™ & (4) 1.2mL prefiled high pressure syringes		
Vista-Cal™ Value Pack	[502601]	\$80.75
(20) 1.2mL prefiled high pressure syringes		

Vista products are available through your preferred dealer

ENDODONTICS

Syringe Warmer

Perfect for heating solutions for WARM delivery of irrigants to surgical site.

- Heats up to 130° F ideal for irrigating solutions
- Cool touch technology housing is comfortable to touch and will not heat other metal surfaces
- Two additional tip recapping slots
- Finger grip side slots
- Runs on standard line voltage

Syringe Warmer	110V	220V	PRICE
6сс	[404301]	[240301]	\$378.75
12cc	[404302]	[240302]	\$378.75
Warmer Bowl Inserts	QTY	SKU	PRICE
	30 / BOX	[404317]	\$41.00



Versa-B<u>rush</u>'

Compare to Ultradent®*: StarBrush®* Rotary Brush

Rotational Utility Brush

Ideal for enhanced cleansing of endodontic access cavities and for removing contaminants.

- For cleaning post holes and coronal access openings
- · Cleansing and removal of smear layer
- · Use with EDTA

QTY	SKU	Price
40 / CUP	[204300]	\$135.15



Citric Acid

Compare to Ultradent®*: Citric Acid

12.5% Solution

Mild chelating agent for removal of smear layer.

- Irrigates, cleans and conditions root canals
- Dissolves calcium hydroxide from canals

	SKU	PRICE
Bulk (30mL) syringe	[503600]	\$22.75
(1) prefilled syringe and (1) syringe docking port		





Vista-Blue™

Methylene Blue

Detects fractures and hard-to-find root canal orifices. Available in 2 sizes for multiple applications of Methylene Blue.

- Unit dose pack
- 10mL bottle kit
- · Individually packaged Unique vial with applicator brush

and (1) 10mL bottle of Methylene Blue

	QTY	SKU	PRICE	
Unit Dose	30 / BOX	[317007]	\$40.45	
Fracture Detection Kit	1 KIT	[503429]	\$37.40	
(24) Micro-Applicator brushes, (6) mixing wells,				



Unit dose delivery Individually packaged unique vial with applicator brush

Micro Aspirators

Compare to Ultradent*: Luer HVE

Designed to provide fast and efficient micro-aspiration.

- · Virtually eliminates the need for paper points
- Pre-tipped style includes 3/4" bendable tips (18, 20 or 25 ga)
- · Autoclavable style can be used with both universal high-or low-speed valve fitments
- Pre-tipped Luer HVE Adapters





COLOR	GAUGE	PKG	SKU	PRICE
	18ga	20 / BAG	[315258]	\$20.30
	20ga	20 / BAG	[315252]	\$20.30
	25ga	20 / BAG	[315255]	\$20.30
Autoclavable :	Style	8 / BOX	[315240]	\$36.50
Luer HVE Unti	pped Adapters	20 /BAG	[315250]	\$16.65
Micro Aspirato	or Evac Set	24 / BOX	[315230]	\$48.30
(1) untiped luer F	HVE adapter			

Fill Station™

Patented

A simple, convenient way to fill syringes!

No mess ... no waste ... less time!

- Simple, convenient way to fill syringes with solutions
- Connect several units to easily fill multiple solutions
- · Syringe-activated valve and low level light indicator
- · Easy to wall mount or use as tabletop modular unit



EZ-Lock™

Patented

A "MUST" for all dental disciplines! The ideal irrigation accessory for subgingival irrigation in restorative, hygiene, periodontics, endodontics and oral surgery.

- · As LOW as \$0.22 per procedure!
- Easily converts any A/W syringe into a subgingival irrigation handpiece when fitted with a luer lock style irrigation tip.
- Small (2") design for better visibility and greater access during irrigation and drying.
- Easy-to-use tips are disposable.
- Eliminates concern for cross-contamination
- EZ-Lock™ is compatible with current Sani-Tip® adapters. Vista adapters are available for most air/water syringes, including A-dec, Dental EZ, Star style, Engle/Marco style, and Press Ring style.



QII	JKU	FRICE.
100 / BAG	[315905]	\$23.40
200 / BAG	[315910]	\$44.50

DDICE.






SKU	PRICE
[315387-A]	\$19.00
[315387-B]	\$19.00
[315387-C]	\$19.00

Precision Irrigator™ Compare to Stropko™* Irrigator

+ **SAVE 47%**

Ideal for use with Blue-Flo™ tips

- For precise control of air and water flow
- · Improves cleaning and drying
- · Small tip design for greater visibility
- · Adapts all luer lock style needle tips
- · Available 2" and 3.5" long
- Autoclavable
- · Can be used with most A/W syringes!
- Essential for micro-endodontic and restorative procedures



EZ-Lock™ Pre-Tipped 3/4" Bendable

(1) Precision Irrigator™ 3.5″ length & (5) Blue-Flo™ tips



COLOR	GAUGE	QTY	SKU	PRICE
	18ga	200 / BAG	[315418]	\$62.00
	19ga	200 / BAG	[315428]	\$62.00
	20ga	200 / BAG	[315429]	\$62.00
	22ga	200 / BAG	[315430]	\$62.00
	25ga	200 / BAG	[315431]	\$62.00




Bleach Safe Valve Bone Trap & Filter Brushes Disposable HVE Screens Endodontic Aspirators Frazier Aspirators High Volume Aspirators High-Speed Suction Valve HV Saliva Ejector Low-Speed Saliva Ejectors Low-Speed Saliva Ejector Valve **Micro Aspirators** Pistol-Grip High-Speed Valve Aspirator **Replacement Heads Saliva Ejector Replacement Filters Straw Aspirator Stylets Surgical Aspirators Surgical Evacuation Tubes Surgical No-Clog**

High-Speed Suction Valve

Autoclavable, quick-disconnect, and flip-style valve

- Autoclavable
- On/off flip action with highly durable metal flip valve
- · Simple to remove and clean flip valve



	SKU	PRICE
Equivalent to: Flip A	[803001]	\$64.65
Equivalent to: Flip D	[803002]	\$64.65

Straw Aspirator

Compare to Pelton & Crane: Tip-A-Dilly®* Aspirator

High-speed suction aspirator

- LENGTH: 4.5"
- 1.5mm opening
- Fits all varieties of straw aspirators
- · Side-venting suction adjust-ment holes regulate aspiration during procedures



Strt	ivvs.	1.5	

PRICE

	Q	51.10	
Straw Aspirator Kit	1 / PACK	[802017]	\$89.50
(1) Aspirator & (200) replacement straws			

Replacement Straws 200 / BOX

[802020] \$12,10

Pistol-Grip High-Speed Valve Aspirator

For Flip A or Flip D style HVE valves

- Autoclavable
- 90° design for reduced hand and arm fatigue caused by conventional aspirator hook-ups

Adapter Style	SKU	PRICE
Standard	[803006]	\$65.30
Bull Frog	[803016]	\$65.30




Surgical Evacuation Tubes

For micro-aspiration use white (1/8")

For macro suction use green (1/4") or blue (3/8")

- Pre-bent to a 30° angle for greater access
- Autoclavable
- •25/bag
- LENGTH: 6"



COLOR	DIAMETER	VENTED	SKU	PRICE
	1/8"	Non-Vented	[403201]	\$12.75
	1/4"	Non-Vented	[403202]	\$12.75
	3/8"	Vented	[403213]	\$12.75
	3/8"	Non-Vented	[403211]	\$12.75

Micro Aspirators

Designed to provide fast and efficient micro-aspiration

- · Virtually eliminates the need for paper points
- Pre-tipped style comes with 3/4" bendable tips (18, 20 or 25 ga)
- · Autoclavable style can be used with both universal

high-or low-speed valve fitments

	COLOR	GAUGE	QTY	SKU	PRICE
Pre-tipped Luer Adapters		18ga	20 / BAG	[315258]	\$20.30
Pre-tipped Luer Adapters		20ga	20 / BAG	[315252]	\$20.30
Pre-tipped Luer Adapters		25ga	20 / BAG	[315255]	\$20.30
Autoclavable Style			8 / BOX	[315240]	\$36.50
Luer HVE Untipped Adapters			20 / BAG	[315250]	\$16.65
Micro Aspirator Evac Set			24 / BOX	[315230]	\$48.30
(1) untiped luer HVE adapter and 40" (extension				








High Volume Aspirators

All Stainless Steel Construction

- Length: 6" / 4"
- Tissue relief slots prevent grabbing of tongue, cheek or soft tissue
- Standard high-volume aspirators
- Bend allows for greater comfort during aspiration



Length: 6"






- Length: 6"
- Tissue relief slots prevent grabbing of tongue, cheek or soft tissue
- Flattened surface to help with retraction




Length: 6"





- · Length: 6"
- Back-cut high-volume aspirator
- Tissue relief slots prevent grabbing of tongue, cheek or soft tissue
- Flattened surface to help with retraction







- Length: 4"
- Short-cut, high-volume aspirator with overall length of 4-1/4".
- Tissue relief slots prevent grabbing of tongue, cheek or soft tissue.
- Flattened surface to help with retraction





OPENING	SAME AS	SKU	PRICE
6 X 4 mm	44FSA	[802013]	\$33.05
8 X5 mm	54FSS	[802014]	\$33.05
9.5 X 6.5 mm	64FSS	[802015]	\$33.05
11.5 X 7.5 mm	74FSSS	[802016]	\$33.05



Disposable HVE Screens

Use with high-volume aspirators to trap debris

- Disposable and cost-efficient
- Adapts to plastic or metal high-volume aspirators
- Helps prevent grasping of soft tissue and stops amalgam

QTY		SKU	PRICE
100 / B	AG	[803011]	\$18.20

Bleach Safe Valve

Bleach-resistant valve



- Bleach-safe on/off flip action
- Longer life cycle and lasting durability

QTY	SKU	PRICE
1 / BAG	[803012]	\$141.45



Endodontic Aspirators

Stainless Steel Construction

Micro suctioning aspirators - ideal for root canal aspiration.



- Finger rest button for exact control of aspirator
 - Includes 1 cleaning stylet and fits HVE cutoff valve
- 1/8" diameter stainless steel
 - Fits HVE cutoff valve
- 1/8" diameter stainless steel.
 - Tapered tip prevents clogging.
 - Fits HVE cutoff valve
 - Bendable tubing for accurate placement.
 - Finger rest button for exact control of aspirator
 - Includes 1 cleaning stylet and fits HVE cut-off valve

OPENING	SAME AS	SKU	PRICE
0.5 mm	01BA	[804001]	\$76.95
1.0 mm	02BA	[804002]	\$76.95
2.0 mm	03BA	[804003]	\$76.95
2.3 mm	05A	[804004]	\$33.05
1.7 mm	04A	[804005]	\$33.05
1.7 mm	04-90A	[804006]	\$33.05
1.7 mm	E04	[804007]	\$33.05
1.7 mm	E04-90	[804008]	\$33.05
1.0 mm	02BAT	[804009]	\$78.55



LENGTH: 6"

804004

804005

804006





804010

LENGTH: 2.75"





Frazier Aspirators

03BAT

Surgical suction aspirator


2.0 mm

Frazier style aspirators fit surgical tubing or Flip A and Flip D high-speed suction valves (#803001 and 803002).

[804010]

\$78.55

- Tapered tip helps prevent clogging.
- · Uses medium size cleaning brush
- Finger suction adjustment hole regulates aspiration during procedures



Low-Speed Saliva Ejector Valve

Autoclavable / Fits standard 3/16" vacuum tubing

- On/off flip action with durable metal flip valve
- Quick disconnect from tubing
- · Simple to remove and clean flip valve

SAME AS	SKU	PRICE			
Flip SE	[803003]	\$61.15			
(1) valve, (3) autoclavable heads,					
and (1) saliva filter					



Flip Style

Replacement Heads

Autoclavable



QTY	SKU	PRICE
3 / BAG	[803004]	\$22.45

Saliva Ejector Replacement Filters

• Fits Vista Saliva Ejector Valve and most other saliva ejectors

QTY	SKU	PRICE
25 / BAG	[803005]	\$6.10
100 / BAG	[803015]	\$12.70



Bone Trap & Filter


For collecting bone fragments, osseous materials and root tips.

- Autoclavable
- Comes with 3 sterile filters
- Fits high-volume cut-off valve or 1/4" tubing



Stylets

- Autoclavable
- · Sturdy and flexible
- Choice of sizes for a variety of aspirator and valve diameters



SIZE	SKU	PRICE
1.0 mm	[900007]	\$7.45
1.5 mm	[900008]	\$7.45
3.0 mm	[900009]	\$7.45

Low-Speed Saliva Ejectors

Use with low speed surgical suction valves (saliva ejectors) and is Autoclavable.



ANGLE	LENGTH	OPENING	SKU	PRICE
30°	4"	1.5 mm	[806001]	\$28.40
30°	4"	2.5 mm	[806002]	\$28.40
30°	4"	3.0 mm	[806003]	\$28.40
90°	4"	1.5 mm	[806004]	\$28.40
90°	4"	2.5 mm	[806005]	\$28.40
90°	4"	3.0 mm	[806006]	\$28.40
30°	6"	1.5 mm	[806007]	\$28.40
30°	6"	2.5 mm	[806008]	\$28.40
30°	6"	3.0 mm	[806009]	\$28.40
90°	6"	1.5 mm	[806010]	\$28.40
90°	6"	2.5 mm	[806011]	\$28.40
90°	6"	3.0 mm	[806012]	\$28.40

HV Saliva Ejector

High-volume suction aspirator

· Adapts to standard high-volume valves






is recommended

Brushes

- Autoclavable
- Choice of sizes for a variety of aspirator and valve diameters
- · Sturdy and flexible

	OPENING	12/PACK	PRICE
Extra Small	3/32"	[900001]	\$37.20
Small	1/8"	[900002]	\$37.20
Medium	1/4"	[900003]	\$37.20
Large	1/2"	[900004]	\$37.20
Extra Large	3/4"	[900005]	\$37.20
Surgical	3/8"	[900006]	\$37.20
Mixed Pack -	2 ea. of above	[900010]	\$37.20



Surgical Aspirators

All Stainless Steel Constuction

• Fits HVE cut-off valve



OPENING	SAME AS	SKU	PRICE
1.5 mm	15P1A	[801001]	\$33.05
2.5 mm	15P2A	[801002]	\$33.05
3.0 mm	15P3A	[801003]	\$33.05
4.0 mm	16RSAA	[801004]	\$33.05
1.5 mm	15P1-90A	[801005]	\$33.05
2.5 mm	15P2-90A	[801006]	\$33.05
3.0 mm	15P3 -90A	[801007]	\$33.05
4.0 mm	16RS-90A	[801008]	\$33.05
1.5 mm	46P1A	[801009]	\$33.05
2.5 mm	46P2A	[801010]	\$33.05
3.0 mm	46P3A	[801011]	\$33.05
2.5 mm	56P2A	[801012]	\$33.05
3.0 mm	56P3A	[801013]	\$33.05
4.0 mm	56P4A	[801014]	\$33.05
1.5 mm	HU1	[801015]	\$63.55
2.5 mm	HU2	[801016]	\$63.55
3.0 mm	HU3	[801017]	\$63.55

1.5mm 2.5mm 3mm 4mm



801006

801007

801008

801009

801010

801011

Length: 6"

- Tapered tip helps reduce clogging in the tube
- Tissue relief slots prevent grabbing of soft tissue (#801004)



- Right angle ideal for socket and endo suction
- Tapered tip helps prevent clogging
- Tissue relief slots prevent grabbing of soft tissue (#801008)



- Made of 1/4" stainless steel
- Large size tube for higher volume suctioning
- Tapered tip helps prevent clogging



Length: 6"

- Made of 5/16" stainless steel
- Large size tube for higher volume suctioning
- Tapered tip helps prevent clogging



Length: 9.25"

 Surgical suction tip with built-in Frazier style vacuum control





3-piece function allows for maximum asepsis

• Fits HVE cut-off valve and uses medium

• Tapered tip to help prevent clogging

control of internal components with built-in stylet.

Surgical No-Clog

size cleaning brush



ANGLE	OPENING	SKU	PRICE
30°	1.5 mm	[801018]	\$115.35
30°	2.5 mm	[801019]	\$115.35
30°	3.0 mm	[801020]	\$115.35





Amdent* Style General Scaling Tips

Amdent* Style Restorative and Condensation Tips

Apical Surgery Tips

AutoSyringe™

Brite-N-Lite™

Cavi-Pro™ Ultrasonic Tips

Co-Pilot™

Composite Applicator™

Endodontic Tips Condensation Tips

Endodontic Tips Restorative Tips

Endodontic Tips Retreatment Tips

EndoUltra™

General Scaling Tips

Periodontal Tips

Pilot™

Step Down Tips

Therma-Flo™

Ultra-Flo™


Ultrasonic Gracey Mini Curette Tips

200X More Powerful the EndoActivator®!



SUPERIOR ACTIVATION* (compared to sonic devices) The EndoUltra™ provides the ultrasonic energy required to effectively produce cavitation

The EndoUltra[™] provides the ultrasonic energy required to effectively produce cavitation resulting in significantly better endodontic procedures.



	SKU	PRICE	
Complete EndoUltra™ Kit	[407702]	\$660.00	
Includes: (1) cordless handpiece (6) 15/02 Activator tips	•	eves,	
	PACK	SKU	PRICE
Activator tips 15/02	(3)	[407706]	\$54.00
Autoclavable Sleeves	(2)	[407709]	\$32.00
Protective Barriers	(100)	[407719]	\$12.00
EndoUltra™ Tip Wrench	(3)	[407711]	\$27.00

COMPARE:	E ndo U ltra [™]	Dentsply®* Tulsa Dental Specialties EndoActivator®
Scientifically Supported Ultrasonic Device	YES	NO
Frequency	ULTRASONIC	SONIC
Multi-use Tips	YES	NO
Cordless	YES	YES

Vista products are available through your preferred dealer

Evidence Supports ULTRASONIC Technology

The next evolution of irrigant activation: ENDO ULTRA

The world's first CORDLESS ultrasonic device!

Research Shows

Conventional hand instrumentation leaves as much as 35% of the canals untouched. Irrigation, although an important part of the canal disinfection and debridement process, cannot in itself be relied on to optimize canal cleanliness.

Proper irrigation protocol should include the use of enhanced endodontic solutions to facilitate removal of microorganisms, tissue remnants, and sentin chips from the root canal. Vista Dental offers a **complete line of enhanced irrigation solutions** that have been tested and proven for superior outcomes.

Irrigants are more effective when they are electro mechanically activated. Acoustic streaming and cavitation have been proven to significantly enhance cleansing of difficult anatomy. Ultrasonic energy not only disrupts biofilm while improving penetration of irrigants and removal of dentinal debris, but it also greatly reduces bacteria levels and improves root sealing.

Ultrasonic energy not only disrupts biofilm while improving penetration of irrigants and removal of dentinal debris, but it also greatly reduces bacteria levels and improves root sealing.

Irrigation Promo Pack

A winning combination for optimal canal cleanliness and improved outcomes!
The only **TRUE** 2-step irrigation protocol
No rinse needed when using Chlor-XTRA™
in conjunction with SmearOFF™ 2-in-1!

Chlor-XTRA™

Enhanced 6% sodium hypochlorite solution

- 2X more dissolution
- 2.5X thinner than standard NaOCL

SmearOFF™

- · Optimizes smear layer removal
- Clears 27% MORE dentin tubules (compared to leading 17% EDTA product)
- · Kills 99.99% bacteria in 10 seconds1
- · Won't form a precipitate when mixed with NaOCI



SKU	PRICE
[317113]	\$79.00

Ultra Efficient. Better cleansing in less time!*

Ultra Effective.

Better penetration of irrigants throughout canal!

Ultra Sonic

Better acoustic streaming with cavitation!

LED light for improved visualization

Autoclavable
15/02 Activator Tips

Contra-angle design for easy access

Removable Autoclavable Sleeve

CALL TOLL FREE 877.418.4782 OR VISIT VISTA-DENTAL.COM





Engineered and manufactured in the U.S.A.

Ultrasonic vs. Sonic Studies Prove:


• That passive needle activation and even low frequency (Sonic) oscillation (160-190Hz) is simply not sufficient in producing acoustic streaming or cavitation within the canal space.

"The positive relationship between streaming velocity and frequency can explain the higher efficiency of PUI [(passive ultrasonic irrigation)] versus sonic irrigation"

Passive ultrasonic irrigation of the root canal: a review of the literature van der Sluis 2007

 That activation by means of ultrasonic frequency (30 kHz) not only generates increased acoustic streaming and cavitation, but results in significantly improved debridement of canal spaces and improved penetration of irrigants into dentinal tubules.

TOOTH CROSS-SECTION Showing Cleaning Efficacy of Canals and Isthmus



Percentages = amount of tissue dissolved

Green = remaining pulp tissue

Study data: ‡Reprinted from Gutarts R, Nusstein J, Reader A, Beck M.
In vivo debridement efficacy of ultrasonic irrigation following hand-rotary instrumentation in human mandibular molars. J Endod. 2005, with permission from Elsevier.



"I feel that EndoUltra™ with lots of irrigation is the only way I got sealer past this 90 degree turn in the distal root."

- Nick Southall, DMD. Mobile, AL






Courtesy of: Dr. Heinrich Dippenaar, DDS. Welkom, South Africa

Chlor-XTRA[™] multiple patents pending









"EndoUltra blends the cordless convenience and the proven adjunctive method of canal cleansing. Once you go Ultra, you won't accept anything less."

Dr. Tung Bui, DDS, FICD, Tuscon, AZ



"EndoUltra offers the most convenient method to deliver passive ultrasonic irrigation. Right weight, portable, ergonomic, efficient, and re-chargeablewhat's not to like?"

- Dr. Joeseph D. Maggio, DDS. Lisle, IL



"I use the EndoUltra unit with NaOCI to disinfect the canals. I am happy with it compared to other units such as EndoActivator®"

- Jamshid Faghih DDS. Bellevue, WA

Smear*OFF*[™] patent pending 2-in-1 mix removes smear layer AND kills 99.9% of bacteria at the same time! **SCAN CODE FOR DATA** To order, see page 33



Ultrasonic Energy Disrupts Biofilm and:

- IMPROVES removal of dentin debris
- IMPROVES penetration of irrigants
- IMPROVES cleanliness
- REDUCES bacteria levels
- IMPROVES accessibility
- IMPROVES root sealing

Ideal for use with Vista engineered endodontic solutions for optimal cleansing and improved results.

EndoUltra™ Disrupts biofilm and:

- IMPROVES cleanliness
- REDUCES bacteria levels
- IMPROVES penetration of irrigants
- Does NOT engage or remove tooth structure
- Resonates down the entire length of the Activator tip
- Runs at 40,000 Hz (vs - sonic versions that run at 160Hz - 190Hz)





Irrigation and Ultrasonic Activation, Together At Last.

Vista Dental Products introduced **UltraFlo™**. This patented device combines the two most vital components of **SUCCESSFUL endodontic treatment**; <u>Irrigation & Ultrasonic Activation!</u>



Advancements in rotary instrumentation have led to shortened contact time of irrigating solutions within the canal system. For this reason, it is more important than ever to properly irrigate. Science has shown that irrigating solutions are more effective upon the introduction of ultrasonic frequencies. This leads to acoustic streaming and cavitation; known to significantly improve canal cleanliness.

UltraFlo[™] is a cordless, compact, battery operated piezo ultrasonic (40kHZ) device that simultaneously irrigates and activates. Providing a constant, consistent flow of endodontic solutions has never been easier. UltraFlo[™] features a quick-connect, 20mL reservoir that is easily filled with the medicament of one's choice. UltraFlo[™] ultrasonic activator is another example of Vista Dental Product's dedication to more safely and effectively achieve endodontic success.



On / Off Button



Autoclavable Sleeve

Flexible, Autoclavable Activator Tip

Ultrasonic Energy Disrupts Biofilm and:

- IMPROVES removal of dentin debris
- IMPROVES penetration of irrigants
- IMPROVES cleanliness
- REDUCES bacteria levels
- IMPROVES accessibility
- IMPROVES root sealing

Ideal for use with Vista engineered endodontic solutions for optimal cleansing and improved results.



8 Tips Included

Piezo Pilot[™] scalers offer easy portability, complementing a light-weight, compact design.

Now with lighted LED handpiece for greater accessibility



Includes (8) ultrasonic tips (1 each) SU1, SU6, HC and PTU, plus 4 tips of your choice.

Pilot[™] Includes:

- 1 Detachable handpiece
- 1 Omni-directional foot pedal
- 1 Quick-disconnect water line
- 3 Autoclavable handpiece sleeves
- 8 Ultrasonic tips
- 8 Torque wrenches

without **LED:** \$1,495

Acteon [402701] EMS [402711]

with **LED:** \$1,695 Acteon [402800]

PIEZO **PILOT**™:

Innovative. High-end performance.

- Push button controls
- Variable irrigation flow settings
- Broad range of power for any endodontic procedure

Versatile design

- · Left or right mounting magnetic handpiece holder
- Detachable handpiece (available in Acteon* or EMS*)

Compact. Comfortable.

- Easy to remove autoclavable silicone sleeves
- Small, lightweight and portable

COMPARE BEST VALUE ¹	VISTA LED Piezo Pilot™	Satelec®* P5 Newtron® LED	Hu-Friedy®* Piezon® 150 LED	Ultradent®* Ultrawave™
Tips Included	8	3	3	3
HandpieceSleeves	3	1	1	1
Price: without LED	\$1,495	\$1,917 ¹	N/A	\$1,799
Price: with LED	\$1,695	\$2,583 ¹	\$1,990¹	N/A

based on competitor's 2015 pricing. *Not a registered trademark of Vista Dental Products.













Vista products are available through your preferred dealer





Piezo Pilot™ scalers offer easy portability, complimenting

a light weight and compact design.

This is a great feature for practices where space is an issue or the unit will be shared.



Includes (10) ultrasonic tips (2 each) SU1 and SU6, (1ea) HC and PTU, plus 4 tips of your choice.

Piezo Co-Pilot™ Co-Pilot™ Includes:

- 1 Piezo Pilot™ ultrasonic scaler
- 1 Detachable handpiece
- 1 Omni-directional foot pedal
- 1 Quick-disconnect water line
- 10 Ultrasonic tips
- 10 Torque wrenches

without LED: \$2,459

Acteon [402901] EMS [402911]

with LED: \$2,599
Acteon [402801]

Piezo Co-Pilot™ Dual Reservoir Scaler

Features (2) 300mL fluid capacity reservoirs - one on each side, with easy-to-see liquid levels



See a video demonstration and an independent review of the Pilot™ at http://cpsmagazine.com/products-unit/piezo-pilot

COMPARE BEST VALUE ¹	VISTA LED Piezo Co-Pilot™	Satelec®* P5 Newtron® XS LED	Hu-Friedy®* Piezon® 250 LED	Ultradent®* Ultrawave™ XS LED
Tips Included	10	3	3	3
Handpiece Sleeves	3	1	1	1
Reservoir Size	600ml	300ml	350ml	300ml
Price: without LED	\$2,459	\$2,554 ¹	N/A	N/A
Price: with LED	\$2,599	\$3,100¹	\$2,890¹	\$2,649

¹based on competitor's 2015 pricing. *Not a registered trademark of Vista Dental Products.













CALL TOLL FREE 877.418.4782 OR VISIT VISTA-DENTAL.com



Cavi-Pro™ Ultrasonic Tips

Ideal for creating small volume cavities for conservative dentistry. For Minimally Invasive Cavity Prep.

- Vista ultrasonic tips are designed to fit all major piezo scaling units including the Piezo Pilot™ Offer direct access to the lesion for greater precision.
- Will not damage adjacent tooth structure or surrounding healthy tissue.
- Minimally invasive for greater patient comfort

Cavi-Pro[™] Tips for Piezo Pilot[™] and...

EMS* STYLE: Mini Piezon*, Scalex 880*, Bonart*, SuperSon EIE*,

Master 400*, Parkel (Piezo)*, Dent America*

Acteon* STYLE: Vista P5*, NSK*, Sybron Endo*, Obtura Spartan*, ADEC*, J.Morita*, Dentsply/Tulsa*, Amdent* Call for tip compatibility with other units.

Price per tip: \$101.45



CP-1 - Half Ball (45 degree right)

For prep, opening, buccal/lingual/occlusal lesions



CP-2 - Half Ball (up)

Use on interproximal lesions and to help preserve adjacent tooth structure



CP-3 - Half Ball (down)

Use on interproximal lesions and to help preserve adjacent tooth structure



CP-P - Curved (rounded conical)

For prep, opening and for buccal/lingual/occlusal lesions



CP-R - Half Ball (45 degrees right)

For interproximal lesions and to help preserve adjacent tooth structure



CP-L - Half Ball (45 degrees left)

Use on interproximal lesions and to help preserve adjacent tooth structure






If your tip is starting to look like this...

Vista has an on-line tip guide to help you maintain the best possible results with your ultrasonic tips.

Endodontic + Condensation tips for Piezo Pilot™ and...

EMS* STYLE: Mini Piezon*, Scalex 880*, Bonart*, SuperSon EIE*, Master 400*, Parkel (Piezo)*, Dent America*

Acteon* STYLE: Vista P5*, NSK*, Sybron Endo*, Obtura Spartan*, ADEC*, J.Morita*, Dentsply/Tulsa*, Amdent* Call for tip compatibility with other units.



Vista products are available through your preferred dealer



Ultrasonic Gracey Mini Curette Tips


Patent Pending

The shorter, more slender blades, plus a longer shank allow these patent pending Gracey Mini Curette tips better access in places that are out of reach with regular sized tips (for instance, with misalignments, areas of tight tissue attachment when furcations are present, or when teeth overlap and are crowded because the bite was never corrected with braces).

- Longer shank than regular curettes: For better access into deeper pockets.
- **Slender blades:** Aid in tactile sensitivity, enhanced deep pocket access and adaptation to root surfaces.
- Shorter blades (average 1.2 mm): For better scaling in NARROW pockets, reduced tissue damage and greater adaptation to the root surface.


Price per tip: \$125.90

BUY 2 or more GET 1 FREE BUY 5 or more GET 3 FREE



For anterior teeth, premolars and molars

Each anatomically specific, uniquely angled Gracey Mini Curette tips adapts EXACTLY to the root surfaces in anterior, molar and pre-molar areas.



Mini Curette Tips for Piezo Pilot™ and...

EMS* STYLE: Mini Piezon*, Scalex 880*, Bonart*, SuperSon EIE*, Master 400*, Parkel (Piezo)*, Dent America*

Acteon* STYLE: Vista P5*, NSK*, Sybron Endo*, Obtura Spartan*, ADEC*, J.Morita*, Dentsply/ Tulsa*, Amdent*

Call for tip compatibility with other units.



HC-1 / HC-2 (anterior teeth)

- Angled for exact adaptation to root surface of anterior teeth
- HC-Mini-1: for mesial/buccal and distal/lingual surfaces
- HC-Mini-2: for distal/buccal and mesial/lingual surfaces



HC-3 / HC-4 (pre-molars)

- Angled for exact adaptation to pre-molar root surface
- HC-Mini-3: for mesial/buccal and distal/lingual surfaces
- HC-Mini-4: for distal/buccal and mesial/lingual surfaces



HC-7 / HC-8 (for molars)

- Angled for exact adaptation to root surface of molars
- HC-Mini-7: for mesial/buccal and distal/lingual surfaces
- HC-Mini-8: for distal/buccal and mesial/lingual surfaces






General Scaling Tips

Our most popular light to heavy-duty tips **Price per tip: \$98.50**

General Scaling Tips for Piezo Pilot™ and...



EMS* STYLE: Mini Piezon*, Scalex 880*, Bonart*, SuperSon EIE*, Master 400*, Parkel (Piezo)*, Dent America* **Acteon* STYLE:** Vista P5*, NSK*, Sybron Endo*, Obtura Spartan*, ADEC*, J.Morita*, Dentsply/Tulsa**, Amdent* **Call for tip compatibility with other units.**





Vista products are available through your preferred dealer




Periodontal Tips

"H" series (hygiene) and "PT" series (perio thin)

Price per tip: \$98.50

Periodontal Tips for Piezo Pilot™ and...







HC - Universal curette design (comes with tip wear marking)

Designed for deep pocket access. Ideal for sub- and supra-gingival scaling.
 Powerful enough for medium to heavy calculus removal. Use edge of tip to quide motion for smooth results.

BUY 2 or more GET 1 FREE



HA - Light to medium sub-gingival calculus (Diamond-coated)

 Use this probe at low power levels to treat anterior teeth or as follow-up to more abrasive sub-gingival tips.



PTU - Sub-gingival pocket scaling (comes with tip wear marking)

Thin tip for sub-gingival scaling; simulates natural curvature of a tooth to provide easy access into pockets with minimal patient discomfort. This universal tip is extremely thin to minimize patient sensitivity and create a smoothly scaled surface.



PTL / PTR - Sub-gingival scaling (PTL - Left Bend / PTR - Right Bend)

Used for lower-left and upper-right molars and pre-molars.
 Also effective for scaling interproximal spaces between molars.



HFL / HFR - Light scaling of difficult areas (Diamond-Coated) (HFL - Left Bend / HFR - Right Bend)

Use on furcations and interproximal sides of molars and pre-molars. Proper use requires non-foot pedal actuated placement, followed by activating foot pedal, scaling, releasing foot pedal and removal of tip. Use at low power settings.













THREAD STYLE	■ HA	■ HC	■ HFL	■ HFR	■ PTL	■ PTR	■ PTU
EMS*:	[502102]	[502103]	[502203]	[502204]	[502201]	[502202]	[502104]
Acteon*:	[401102]	[401103]	[401203]	[401204]	[401201]	[401202]	[401104]
Amdent*:	[403102]	[403103]	[403203]	[403204]	[403501]	[403502]	[403104]
Use for:	Light to medium Supra-Gingival Calculus	Universal curette design for sub and Supra-Gingival Calculus	Left Bend light scaling of hard-to-reach areas (furcations and interproximal)	Right Bend Sub-gingival pocket scaling of molars and pre-molars	Left Bend Sub-gingival pocket scaling of molars and pre-molars	Right Bend Sub-gingival pocket scaling of molars and pre-molars	Sub-gingival pocket scaling
	\$98.50	\$98.50	\$98.50	\$98.50	\$98.50	\$98.50	\$98.50



Amdent* Style

General Scaling Tips Price per tip: \$117.75



#33 / #33R / #33L - Probe Tip

For supra-gingival scaling of overhangs and lingual, buccal and interdental scaling.





#37 / #37R / #37L - Thin Tip For fine-supra-gingival scaling, spot removal and sub-gingival

caculus removal.

Amdent* Style

Restorative and Condensation Tips

Price per tip: \$117.75 - \$172.55







Inlay and Onlay
Condensation tips
for inlay and onlay.

Endodontic Tips


THREAD STYLE

Amdent*:

Type of Tip:


PRICE:

Condensation Tips Price per tip: \$98.50





BUY 2 or more GET 1 FREE






Vista products are available through your preferred dealer



Endodontic Tips

Restorative Tips Price per tip: \$101.45

BUY 2 or more GET 1 FREE **BUY 5 or more GET 3 FREE**





Versatile tips are ideal for use in a wet or dry field.

Applications: refining axial walls, caries removal, cleaning pulp chamber, sanding dentin, finishing access cavity, and sanding and finishing restoratives.

-D=Diamond Coated



Endodontic Tips

Retreatment Tips

Price per tip: \$98.50 - \$101.45

BUY 2 or more GET 1 FRE BUY 5 or more GET 3 FREE

> UCT (cavitation tip) 10mm of diamond coating for increased energy transfer. Increases effectiveness of medicaments used during endodontic

PRT (post removal tip) Ideal for loosening of

procedures.

irrigation. Irrigation is recommended as a cooling agent when removing crowns

and posts.



GE20 / GE20D (20mm tip to bend)

Versatile - ideal for use in either wet or dry fields.

Applications: refining axial walls, caries removal, cleaning pulp chamber, sanding dentin, finishing access cavity, and sanding and finishing restoratives.

D= Diamond Coated



GE40 / GE40D (40mm tip to bend)

Versatile - ideal for use in either wet or dry fields.

Applications: refining axial walls, caries removal, cleaning pulp chamber, sanding dentin, finishing access cavity, and sanding and finishing restoratives.

D= Diamond Coated







GE40D

\$101.45



GE40

\$98.50



[502623]

[401623]

None

\$98.50

UCT

[502640]

[401640]

[403640]

Diamond

\$101.45

THREAD STYLE GE20D EMS*: Acteon*: posts with ultrasonic Amdent*: Type of Tip: PRICE: \$101.45

[502301] [502401] [502302] [502402] [401301] [401401] [401302] [401402] [403301] [403401] [403302] [403402] Diamond None Diamond None

\$98.50

CALL TOLL FREE 877.418.4782 OR VISIT VISTA-DENTAL.com



Apical Surgery Tips

Price per tip: \$98.50 - \$101.45

Apical Surgery Tips for Piezo Pilot™ and...

EMS* STYLE: Mini Piezon*, Scalex 880*, Bonart*, SuperSon EIE*, Master 400*, Parkel (Piezo)*, Dent America* Acteon* STYLE: Vista P5*, NSK*, Sybron Endo*, Obtura Spartan*, ADEC*, J.Morita*, Dentsply/Tulsa®*, Amdent* Call for tip compatibility with other units.

BUY 2 or more GET 1 FREE BUY 5 or more GET 3 FREE .

		•			T	T
THREAD STYLE	■ S70D	■ S70	■ S90D	■ 590	SAD	■ SA
EMS*:	[502306]	[502406]	[502303]	[502403]	[502626]	[502636]
Acteon*:	[401306]	[401406]	[401303]	[401403]	[401626]	[401636]
Amdent*:	[403306]	[403406]	[403303]	[403403]		
Type of Tip:	Diamond Coated	Plain Tip	Diamond Coated	Plain Tip	Diamond Coated	Plain Tip
	70° angled tips generally	70° angled tips generally	90° angled tips generally	90° angled tips generally	Universal design generally	Universal design generally
	used on anterior teeth	used on anterior teeth				
PRICE:	\$101.45	\$98.50	\$101.45	\$98.50	\$101.45	\$98.50

THREAD STYLE	SPLD	SPL SPL	SPRD	SPR	■ SMLD	■ SML
EMS*:	[502304]	[502404]	[502305]	[502405]	[502621]	[502634]
Acteon*:	[401304]	[401404]	[401305]	[401405]	[401621]	[401634]
Type of Tip:	Diamond Coated	Plain Tip	Diamond Coated	Plain Tip	Diamond Coated	Plain Tip
	Left angle - generally	Left angle - generally	Right angle - generally	Right angle - generally	Left angle - generally	Left angle - generally
	used on posterior teeth					
PRICE:	\$101.45	\$98.50	\$101.45	\$98.50	\$101.45	\$98.50

THREAD STYLE	■ SMRD	SMR
EMS*:	[502622]	[502635]
Acteon*:	[401622]	[401635]
Type of Tip:	Diamond Coated	Plain Tip
	Right angle - generally	Right angle - generally
	used on posterior teeth	used on posterior teeth
PRICE:	\$101.45	\$98.50



■ S70 / S70D 70° degree angled surgery tip

Generally used on anterior teeth.

The S70D is diamond-coated to offer

greater durability and provide more abrasive

power.



SA / SAD - Universal design surgery tip Generally used on anterior teeth.
SAD is diamond-coated to offer greater durability and provides more abrasive power.



S90 / S90D - 90° degree angled surgery tip

Generally used on anterior teeth.

The S90D is diamond-coated to offer

greater durability and provide more

abrasive power.



■ SML / SMLD - Left angle surgery tip

SMR / SMRD - Right Angle surgery tip

Generally used on molars and

pre-molars or on posterior teeth.


SMLD and SMRD are diamond-coated to offer greater durability and provide more abrasive

power.



SPL / SPLD - Left angle surgery tip
SPR / SPRD - Right angle surgery tip
Generally used on molars and
pre-molars or on posterior teeth.

SPLD and SPRD are diamond-coated to offer
greater durability and provide more abrasive
power.



Therma-Flo®

Revolutionizing Composite Delivery Science Behind the Heat

The heating technology of Therma-Flo™ significantly increases the flow characteristics of highly filled composites while helping to increase polymerization and improve depth of cure. Heat lowers the viscosity of composite material, allowing it to better flow and adapt to cavity walls. As a result, voids are reduced, which then reduces the risk of secondary caries. Heating composite material also improves its physical and handling properties without modifying color or stability.

FEATURES

Benefits of HEAT:

- · Creates fewer voids, with greater depth of cure and shorter curing time
- · There is no adverse effect to pulp tissue or to composite material
- · Easier to place and manipulate composite material



"Incremental placement of composite material becomes extremely accurate with Therma-Flo™. When I use the Composite Warmer with either 2X or 4X Step Down tips, I can easily place precise increments and sculpt the composite into beautiful restorations that will last the patient for years and years."

Randy Shoup, DDS

Heats at tip to maintain optimal delivery temperature	4	


Depth of Cure	Fill restorations faster by increasing the depth of cure and reducing curing time.	⁴ Burtscher 2005, ² Munoz 2008
Flowability	When heated, highly filled composites flow up to 10X better.	³Lucey 2010
Increased Polymerization	Significantly higher monomer conversion values. Dramatic increase in polymerization rates.	⁵ Trujillo 2003, ⁶ Daronch 2005
Decreased Voids	Heating significantly reduces the chance of secondary caries.	⁹ Choudhary 2011, ¹⁰ Wagner 2008, ⁸ Froes-Salgado 2010
Color + Stability	Heating will NOT modify color or stability properties of composite material.	¹Mundim 2011
Micro-Hardness	Heating results in shorter curing times and enhances subsequent surface hardness.	² Munoz 2008, ³ Lucey 2010
Viscosity -vs- Temperature	Heating makes it easier to place material and results in better adaptation to cavity walls.	³Lucey 2010
Heating Safety	Heating will NOT damage pulp tissue or cause discomfort. No other safety concerns.	⁷ Daronch 2007 ¹¹ Rueggeberg 2010


BENEFITS OF USING HEAT

PROOF









Vista products are available through your preferred dealer

Therma-Flo®

Revolutionizing Composite Delivery **Composite Warming Kit**

Increase the flowability of highly filled composites without compromising tensile strength, wear-resistance or aesthetics


- · Base unit warms multiple capsules Includes:
- (1) Warming Unit
- (2) Composite Syringes
- (20) 2X Step Down Tips
- (20) 4X Step Down Tips

VOLTAGE	SKU	PRICE
110V	[404350]	\$595.00
220V	[404355]	\$595.00



Step Down Tips

Extended length: deeper access and precise placement A completely new delivery option enabling doctors to use much less composite with little manipulation required to place the material.



- · Narrow profile for optimal visualization
- · Easy extrusion of composite material

SIZE	QTY	SKU	PRICE
2X	40 / CUP	[320816]	\$30.90
4X	40 / CUP	[320817]	\$30.90

Use your favorite composites!

Therma-Flo™ Applicator™, Warming Kit and Step Down tips are compatible with most manufacturer's capsules.



Composite Applicator™

Patented

Motorized extrusion of heated composites

- · Heats at the point of delivery to maintain optimal flowability throughout procedure
- Improves composite handling properties
- · Lightweight ergonomic design reduces hand fatigue
- Cordless and compact

Includes:

- (20) 2X Step Down Tips
- (20) 4X Step Down Tips
- (1) Rechargeable handpiece
- (1) Warming base

QTY	SKU	PRICE
1	[407821]	\$995.00





Heating up **Composite Delivery**

Vista Dental has revolutionized composite delivery with their new Therma-Flo™ line of products, uniquely engineered to utilize heat for optimal performance of any preferred composite material. Vista's Composite Applicator™ and Composite Warming Kit allow you to use highly filled composites in a wider variety of applications and also provide the versatility to enhance virtually all restorations.

Designed to enhance the performance of Therma-Flo™ products, Vista's Step Down tips offer extended length for deeper access and greater precision in composite placement.

Uniquely designed to adapt to most composite capsules, Step Down tips represent a major step forward in the control of composite delivery. placement.

Uniquely designed to adapt to most composite capsules, Step Down tips represent a major step forward in the control of composite delivery.



SCAN TO SEE HOW IT WORKS!



LIGHT UP YOUR PRACTICE WITH THE NEW BRITE-N-LITETM ORAL ILLUMINATION SYSTEM

Oral Illumination System

Vista Dental™ introduces Brite-N-Lite™, a powerful oral illumination device. Brite-N-Lite™ offers direct illumination of the working area during any procedure. This clip-on LED device is autoclavable and easy to use with virtually any aspirator. The cool white LED light optimizes visibility within the oral cavity.

Other illumination systems on the market require cumbersome battery packs, cords, and expensive adaptors. Brite-N-Lite™ is a compact, convenient illumination system that adapts to the instruments that you currently own. The small, hands-free device securely mounts to your standard evacuation tubing. Brite-N-Lite™ features a long lasting, rechargeable battery for hours of procedural use. The intuitive design of Brite-N-Lite™ makes it the most practical, easy to use oral illumination system on the market.

FREE with every Brite-N-Lite™ comes with one metal aspirator of your choice!

Brite-N-Lite™ LED Kit:

[508001] **\$525.00**

Brite-N-Lite™ LED Light Guide:

[508010] **\$175.00**

Grommets - 5mm:

[508013] **\$25.00**

Grommets - 6.5mm:

[508012] **\$25.00**




Brite-N-Lite™ LED KIT

The illumination kit includes the following components plus one aspirator of choice *Includes*:

- Battery Pack and Charger
- (2) LED Light Guide
- Choice of surgical aspirator
- (12) grommets to fit 5mm and 6.5mm aspirators
- Hose Strap
- Velcro Straps



Visit us at www.vista-dental.com



Air / Water Syringe Tips
Air / Water Syringes
Aspirating Syringes
Bendable Brush Applicators
Bite Wing Tabs
Cheek Retractors
Composite Brushes & Brush Handles
Cotton-Tipped Applicators
Disposable HVE Screens
Lab Knives
Long Shank Brushes
Metal Mirror Stems
Regulated Air / Water Syringes
Spatulas
Syringe Tip Covers

GENERAL MERCHANDISE

Mixing Wells

White, disposable, and able to use with all materials.

- Dual-well style
- Convenient




QTY	SKU	PRICE
100 / BAG	[503414]	\$13.50

Disposable Mouth Mirrors

For office examinations and air abrasion.

Made with High Quality Plated Glass.

- Highly reflective ultra-thin PLATED GLASS surface enhances imaging
- Ideal for high-volume where interim sterilization is not practical
- Light-weight with ergonomic handle
- Can also function as a cheek retractor



QTY	SKU	PRICE
60 / BOX	[307001]	\$27.75

Cheek Retractors

- Disposable
- Retracts cheeks and depresses tongue
- Allows for unrestricted access to the mouth



QTY	SKU	PRICE
10 / PACK	[205215]	\$15.30

Metal Mirror Stems

Fits simple stem or cone socket mirrors.

• Resists Corrosion



SIZE	SKU	PRICE
Round: Simple Stem	[103001]	\$2.25
Cone Socket	[103002]	\$2.25
Octagonal: SImple Stem	[104001]	\$2.25
Octagonal: Cone Stem	[104002]	\$2.25

Spatulas

Wood Handled

Ideal for mixing and spreading

- Durable stainless steel blades
- Solid rosewood handles
- Premium quality





SPATULAS	SKU	PRICE	#8R
#2R	[105020]	\$14.80	πOI
#3R	[105030]	\$14.80	
#4R	[105040]	\$14.80	#10R
#8R	[105080]	\$14.80	
#10R	[105010]	\$14.80	#11
#11R	[105011]	\$14.80	

Lab Knives

For all dental environments - Wood handled

- Solid rosewood or green enamel wood handles
- · Stainless steel, razor-sharp knife blade
- Quality construction





Composite Brushes & Brush Handles

Universal style for application of etchants, sealants, adhesives, and hemostatic agents.

- Tapered tip allows smooth, even application of all liquid types
- DuPont®* nylon bristle tips fit all leading brush handles
- Will not slip or dislodge
- Handle accepts all composite brush styles (including leading brush tips)
- Manufactured in rigid plastic for REPEATED AUTOCLAVING
- Added reach makes application easier

SIZE	QTY	SKU	PRICE
Fine Bristle	100 / BAG	[302001]	\$7.80
Regular Bristle	100 / BAG	[302002]	\$7.80
Handles	2 / BOX	[302004]	\$7.80




Bendable Brush Applicators

For applying of etchants, resins and bonding agents.

- BENDS TO ANY angle and stays in place
- Durable brush with taper-cut bristles
- · Adapts easily for lingual or interproximal use

BENDABLE BRUSH APPLICATORS	QTY	SKU	PRICE
Fine Bristle	144 / BOX	[204001]	\$24.60
Regular Bristle	144 / BOX	[204002]	\$24.60
Combo (half regular & half fine)	144 / BOX	[204005]	\$24.60
Fine Bristle	432/ BOX	[204003]	\$67.35
Regular Bristle	432/BOX	[204004]	\$67.35
Combo (half regular & half fine)	432/ BOX	[204006]	\$67.35

NEW! Bendable Handle Design



Shown with regular bristles.

Cotton-Tipped Applicators

For delivery of all dental liquids and resins.

- Non-sterile virgin cotton
- Tightly woven for lint-free delivery



SIZE	1000/ BOX	PRICE
3 in	[202003]	\$8.80
6 in	[202006]	\$8.80

Long Shank Brushes

For facial, buccal or interproximal applications.

- Black nylon bristles are compatible with ALL dental liquids and resins
- Easy-grip 2-inch handle length

QTY	SKU	PRICE
1000 / BOX	[302003]	\$6.20



Bite Wing Tabs

Identify / pre-tab x-rays

- Sticks indefinitely to all types of bite wing x-ray surfaces
- Unique dispenser box peels off each tab automatically

QTY	SKU	PRICE
500 / Box	[308001]	\$6.99



Disposable HVE Screens

Use with high-volume aspirators to trap debris.

- Disposable and cost-efficient
- · Adapts to plastic or metal high-volume aspirators
- Helps prevent grasping of soft tissue and stops amalgam

_	QTY	SKU	PRICE
	100 / BAG	[803011]	\$18.20



GENERAL MERCHANDISE

Air / Water Syringes

Vista offers premium A/W syringes at non-premium prices that other manufacturers can't beat!

- Autoclavable
- Available in one- and two-button styles
- Includes Fast-Change tip



		1001
AIR / WATER SYRINGES	SKU	PRICE
One-Button (syringe only)	[315390]	\$69.65
(1) One-Button A/W Syringe & (1) Fast Change Tip		
Two-Button (syringe only)		
(1) Two-Button A/W Syringe & (1) Fast Change Tip	[315350]	\$69.65

Syringe Tip Covers

For infection control with metal syringe tips.

- Autoclavable
- Requires no conversion adapters
- · Will not impede air or water flow
- Slides over 3-way metal syringes
- Hydrophobic filter adds
- "suck-back" protection
- Meets OSHA and CDC Standards

QTY	SKU	PRICE
100 / BAG	[301001]	\$17.75
250 / BOX	[301002]	\$39.60



Regulated Air / Water Syringes

With built-in flow regulator

- Autoclavable
- Syringe features a fixed orifice flow control
- Regulated flow on both one and two button syringes
- Ideal for use with EZ-Lock™ and Precision Irrigator



REGULATED A / W SYRINGES	SKU	PRICE
Regulated One-Button Syringe Kit	[315384]	\$109.20
(1) Regulated One-Button A / W Syringe &		
(1) Fast Change Tip and 7' Tubing		
Regulated Two-Button Syringe Kit	[315382]	\$109.20
(1) Regulated Two-Button A / W Syringe &		
(1) Fast Change Tip and 7' Tubing		
EZ-Lock Regulated One-Button Syringe Kit	[315380]	\$178.85
(1) Regulated One-Button A / W Syringe Kit,		
(1) EZ-Lock Adapter, and (200) EZ-Lock™ Tips		
EZ-Lock Regulated Two-Button Syringe Kit	[315378]	\$178.85
(1) Regulated One-Button A / W Syringe Kit,		
(1) EZ-Lock Adapter, and (200) EZ-Lock™ Tips		

Air / Water Syringe Tips

Quick-change

- Autoclavable
- Durable, stainless steel construction



QTY	SKU	PRICE
1 / BAG	[315275]	\$5.65
5 / BAG	[315375]	\$25.05
50 / BAG	[315385]	\$213.15

Aspirating Syringes

Ideal for use with 1.8cc cartridges.

Perfect for warm delivery of anesthetic when used with Vista's 6cc or 12cc Syringe Warmers.

(See page 36 for details on Syringe Warmers.)

- Chrome construction
- Autoclavable
- Durable withstands repeated use

ASPIRATING SYRINGES	SKU	PRICE
Astra Style	[107001]	\$21.25
Cook-Waite Style	[107002]	\$21.25
REPLACEMENT HUBS		
Astra Style	[304001]	\$2.05
Cook-Waite Style	[304002]	\$2.05



MOVE BEYOND THE LIGHT.



A NEW VISION IN THE FIGHT
AGAINST ORAL CANCER HAS ARRIVED.

TOGETHER WE WILL SAVE LIVES.

Find out what 9 years of scientific research has done for discovering oral cancer.

ESPYOC.COM
Visual Solutions LLC.




INDEX		EndoUltra™	
1" Appli-Vac™ Bendable Tips	9, 35	EZ-Lock™	37
1" Secure-Lock™ Appli-Vac™ Bendable Tips		EZ-Lock™ Adapter	37
1" Secure-Lock™ Vista-Probe™ Bendable Tips		EZ-Lock™ Pre-Tipped 3/4" Bendable	37
1"Vista-Probe™ Bendable Tips		Fill Station™	37
1/2" Bendable Applicator Tips		Flexi-Glide™ Utility Tip	8,34
3/4" Appli-Vac™ Bendable Tips		Fluorescent™ Dr. Kit	21
3/4″ Secure-Lock™ Appli-Vac™ Bendable Tips		Fluorescent™ EXpress White™	22
Air / Water Syringe Tips		Fluorescent™ Refill	21
Air / Water Syringes		Frazier Aspirators	41
Amdent* Tips		General Scaling Tips	55
General Scaling Tips	57	High Pressure Syringes	11
Restorative and Condensation Tips		High Volume Aspirators	40
Apical Surgery Tips		High-Speed Suction Valve	
Appli-Vac™ Pre-tipped Syringes		HV Saliva Ejector	
Aspirating Syringes		HVE Tongue Cleaners	
AutoSyringe™		Impression Syringe	
Bendable Brush Applicators		Impression Trays	
Best-Etch™		In-Office Whitening Liquid Dam Kit	
Bite Wing Tabs		Interproximal Brushes	
Bleach Safe Valve		Irrigation Promo Pack	
Blue-Flo™ Tips		Lab Knives	
Bone Trap & Filter		Liquid Dam Refill Kits	
Brushes		Long Shank Brushes	
Brite-N-Lite™		Low-Speed Saliva Ejectors	
CanalClean™ AA		Low-Speed Saliva Ejector Valve	
CanalClean™ CA		Luer Lock Syringes	
Caries Indicator		Metal Mirror Stems	
Cavi-Pro™ Ultrasonic Tips		Micro-Applicator Brushes	
Cheek Retractors		Micro Aspirators	
Chlor-XTRA™		Micro Flowthru® Brush Tips	
CHX Restorative		Micro-Evac™ Tip	
CHX-Plus TM		Micro-Flo™ Tips	
Citric Acid		Micro-Restorative Best-Etch™	
Colored Luer Lock Syringes		Mixing Wells	
Composite Applicator		NiTi SuperFlex™ Tips	
Composite Brushes & Brush Handles		Periodontal Tips	
Cotton-Tipped Applicators		Perioflex™ Irrigator Bottle	
Disposable HVE Screens		Piezo Co-Pilot™	
Disposable Mouth Mirrors		Piezo Pilot™	=
Disposable MTA Carriers		Pistol-Grip High-Speed Valve Aspirator	-
Dual Arch Fluoride Trays		Pocket-Probe™	
EDTA		Porcelain Etch	-
Endodontic Aspirators		Pre-Bent Tips	7
Endodontic 7/5pirators	1 1	Precision Irrigator™	
Condensation Tips	57	EZ-Lock™ Pre-Tipped 3/4″ Bendable	
Restorative Tips		Precision Irrigator™	
Retreatment Tips		Quick-Stat™ FREE	
EndoSafe™		Quick-Stat™ FS	

Quick-Stat'™ Liquid	
Quick-Stat™ TM Bulk Basic Kit	13
Quick-Stat™TM Bulk Complete Kit	13
Regulated Air / Water Syringes	67
Replacement Heads	
Straws Aspirator	39
Safe-T™ Syringes	11
Saliva Ejector Replacement Filters	42
Secure-Lock™ Appli-Vac™	35
Secure-Lock™ Vista-Probe™	
Single Arch Fluoride Trays	18
SmearOFF™	
Sodium Hypochlorite	33
Sol-u-Flo™ Applicator Tips	
Spatulas	65
Spira-Flo™ Brush Tips	8
Stat-Flo™ Tips	8, 14
Stat-Pak™	
Stat-Pak™ +E	
Step Down Tips	62
Straw Aspirator	
Stylets	42
Surgical Aspirators	43
Surgical Evacuation Tubes	39
Surgical No-Clog	
Syringe Caps	11
Syringe Tip Covers	67
Syringe Warmer	
Therma-Flo™	
Tongue Cleaners	18
Ultrasonic Gracey Mini Curette Tips	54
Ultra-Flo™	49-50
US-12 Utility Syringe	11
Versa-Brush [™]	
Visco-Tip™	
Vista Spring™ Tooth Brushes	
Vista-Blue™	
Vista-Breeze™ Toothbrushes	17
Vista-Cal™	35
Vista-Pak™ Packers	
Vista-Probe™ Pre-tipped Syringes	
Vista-Probe™ Tips	
Volu-Flo™ Tips	8
Warmer Bowl Inserts	



2200 South St. Racine, WI, 53404



Call toll free 1.877.418.4782 or go online to www.vista-dental.com